

# Product Catalogue

WINSIX PRODUCTS



<b>COMPANY</b>	
Customer Support and After Sales Service .....	pag. 4
WINSIX Implant System Registered Trademarks .....	pag. 5
Company - Quality - Certification - Patents - Insurance	
Qualified Training - Patient Information .....	pag. 6 - 7
Summary Table: WINSIX Implant indications .....	pag. 8
<b>WINSIX IMPLANT SYSTEM</b> .....	pag. 9
WINSIX Implants: diameters and lengths .....	pag. 10
KAPPA line .....	pag. 10
Torque Type Line .....	pag. 11
WINSIX Micro Rough Surface .....	pag. 12
Implant Collar Types and Heights .....	pag. 13
WINSIX Implant Connections Prospectus .....	pag. 14
Stable fit Implant - Abutment .....	pag. 15
Implant packaging - Opening of packaging - Extraction of the implant .....	pag. 16
<b>KAPPA IMPLANT LINE</b> .....	pag. 17
Features .....	pag. 18 - 19
K Implant .....	pag. 20 - 21
KT Implant .....	pag. 22 - 24
KT MACHINED Implant .....	pag. 25
KE Implant .....	pag. 26 - 27
Soft tissue management .....	pag. 28 - 29
<b>TORQUE TYPE IMPLANT LINE</b> .....	pag. 31
Features .....	pag. 32 - 33
TTi IMPLANT - internal connection .....	pag. 34 - 35
TTx IMPLANT - external connection .....	pag. 36 - 37
TTc IMPLANT - conical connection .....	pag. 38 - 39
TTcm IMPLANT- morse taper connection .....	pag. 40 - 41
<b>SURGERY</b> .....	pag. 43
WINSIX SURGICAL KIT .....	pag. 44 - 45
Equipment available for the Surgical Kit .....	pag. 46 - 47
Summary table for the use of Kappa Line Drills .....	pag. 48
Summary table for the use of Torque Type Line Drills .....	pag. 49
Drills.....	pag. 50 - 51
Bone mill .....	pag. 52
Implant Drill Stops .....	pag. 53
Intraoperative Instruments .....	pag. 54
Intraoperative Accessories .....	pag. 55
Surgical and Prosthetic Instruments .....	pag. 56 - 58
<b>PROSTHETICS - Internal CONNECTION KAPPA LINE - TORQUE TYPE LINE</b> .....	pag. 59
Healing caps .....	pag. 60
Implant closure caps .....	pag. 61
Abutment Impression Copings .....	pag. 62 - 63
Analogues .....	pag. 63
PROSTHETIC ACCESSORIES OVERVIEW .....	pag. 64 - 65
Screw-retained prosthesis .....	pag. 66
Prosthetic Abutments .....	pag. 67 - 74

<b>PROSTHESES: WINZIR AND PRO ABUTMENTS FOR ZIRCONIA PROSTHESIS</b> .....	pag. 75
Prosthetic Abutments .....	pag. 76
<b>PROSTHESES - OVERDENTURE COMPONENTS - Internal CONNECTION</b> .....	pag. 77
LOCATOR Attachments .....	pag. 78 - 79
Ball attachments .....	pag. 80
<b>MULTI-UNIT SCREW-RETAINED PROSTHESIS</b>	
Flat Shift Line - Extreme Abutment EA Line .....	pag. 81
FLAT SHIFT LINE .....	pag. 82 - 89
EA LINE - Extreme Abutment .....	pag. 90 - 97
WINSIX internal hex screws .....	pag. 98
<b>PROSTHESES -TTx TORQUE TYPE EXTERNAL CONNECTION</b> .....	pag. 99
Healing caps .....	pag. 100
Implant closure caps .....	pag. 100
Abutment Impression Copings .....	pag. 101
Analogues .....	pag. 101
Prosthetic Abutments .....	pag.102 - 104
<b>EAX LINE EXTREME ABUTMENT MULTI-UNIT SCREW-RETAINED PROSTHESIS</b> ..	pag.105 - 110
WINSIX external hex screws .....	pag. 111
<b>PROSTHESES - OVERDENTURE COMPONENTS - EXTERNAL CONNECTION</b> ..	pag. 113
LOCATOR Attachments .....	pag. 114 -115
<b>PROSTHESIS - CONICAL CONNECTION</b> .....	pag. 117
Healing caps .....	pag. 118
Implant closure caps .....	pag. 118
Abutment Impression Copings .....	pag. 119
Analogues .....	pag. 119
Prosthetic Abutments .....	pag.120 - 122
Multifunctional Abutments .....	pag.123 - 124
Conical connection transverse screws .....	pag. 124
<b>CAB CLIP ABUTMENT BAR DEVICE</b>	
CAB .....	pag.125 - 130
<b>WINSIX DIGITAL</b> .....	pag. 131
WinScan Abutments .....	pag.132 - 133
<b>WINSIX PROSTHETIC SOLUTIONS</b> .....	pag.134 - 135
<b>MARKETING TOOLS</b> .....	pag. 136
<b>OTHER PRODUCTS</b> .....	pag. 137
WINSIX Flash Implants .....	pag. 138
Osteotomes .....	pag. 139
Safety Zone Kit .....	pag. 139
<b>MATERIALS</b> .....	pag. 140
<b>STERILISATION PROTOCOL</b> .....	pag. 141
<b>MAIN BIBLIOGRAPHY</b> .....	pag.142 - 143



## Customer Support and After-Sales Service

BIOSAFIN pays particular attention to the Technical and Commercial Training of its staff in the area. Our dedicated Product Specialists and Company after-sales service, are available to our customers for technical information, updates and whatever is needed in relation to the marketed products.

**WINSIX Implant System**  
registered in 58 Countries around the world



- |                                 |  |  |
|---------------------------------|--|--|
| 1995 WINSIX®                    | 2010 Teeth Just On 6®                  |  |
| 1998 Free Tense System®         | 2010 Torque Type®                      |  |
| 2001 Bioactive Covering®        | 2011 Clip Abutment Bar CAB®            |  |
| 2001 Free Lock®                 | 2012 WINClinic®                        |  |
| 2007 Full Contact Covering FCC® | 2013 Double Conical Connection DCC®    |  |
| 2009 Flat Shift System®         | 2013 WINPeek® Abutment                 |  |
| 2009 Micro Rough Surface MRS®   | 2016 Extreme Abutment Multifunctional® |  |
| 2010 Extreme Abutment®          | 2016 Linea KAPPA®                      |  |
| 2010 Teeth Just On 4®           |  |  |

- 2009 BioBone®
- 2009 Easy Surgery®
- 2015 Easy Weld®
- 2015 Easy Light®
- 2015 Easy Physio®
- 2017 BiAligner®

Listed trademarks are cited in this edition of the WINSIX Product Catalogue, the omission of the ® is intended to be understood.

WINSIX, Free Tense System, Bioactive Covering, Free Lock, Full Contact Covering FCC, Flat Shift System, Micro Rough Surface MRS, Extreme Abutment, Teeth Just On 4, Teeth Just On 6, Torque Type, Clip Abutment Bar CAB, WINClinic, Double Conical Connection DCC, WINPeek Abutment, Extreme Abutment Multifunctional, Linea KAPPA, BiAligner.

# Specialist Company

**BIOSAFIN** specialises in the production and sale of devices and instruments for Implantology and Oral Surgery, organically combined in the WINSIX Implant System and BIOSAFIN ORAL SURGERY Line.

Research and Development is at the heart of our business, aimed at implementing product performance through the development of innovative technological and scientific content, and constantly improving its Quality. The evolutionary path of the WINSIX System finds scientific evidence in the extensive Bibliography, Studies and Publications available.

## Product quality scientifically and clinically tested.

For over 22 years, the WINSIX Implant System has been used with satisfaction by dental centres and professional studios with varying requirements according to type of consumer and workflow.

The certified quality of the products - which all undergo a strict 1:1 production check - and the proposed solutions, always relevant in meeting the needs of patients, provide the dental surgeon with maximum operational peace of mind.

The solid scientific background upon which devices are developed allows the combination of technological innovation and compatibility, avoiding demanding changes in operating times or expensive equipment replacement for the dental Team.

PRODUCT CERTIFICATIONS:



Certification for the marketing of WINSIX PRODUCTS in USA

## Company certification

**BIOSAFIN** is a Certified Company: **UNI EN ISO 9001: 2015**, which certifies the entire 360° work process, guaranteeing compliance with the quality standards considered to be optimum for the protection of the Product User - the Professional - and the end user - the Patient.

**UNI CEI EN ISO 13485: 2016** specifically relevant to the Quality of Medical Devices.

The quality standards imposed by the Certification are periodically reviewed and re-evaluated.



## Patents: CAB Clip Abutment Bar



International and European Patent  
PCT/EP2011/072448  
EP Patent no. 11425032.7



The patent issued for the CAB by the US authority, confirms its originality and innovative technological content.



The patent issued for the CAB - Clip Abutment Bar by the US authority, confirm its originality and innovative technological content.

## Insurance policies

The WINSIX Implant System devices benefit from TPLI (Third Party Liability Insurance) and CLI (Civil Liability Insurance).

## Qualified Training

BIOSAFIN offers a professional Updating and Training programme in line with corporate quality standards for the correct application of implant devices. The proposed teaching Formats are the result of experience and careful evaluation, differentiated according to the needs of professionals and developed in collaboration with experienced teachers.

REFERENCE CENTRE FOR TRAINING:  
WINSIX Pilot Centre - Vita-Salute San Raffaele University of Milan

Please visit our Training page at [www.biosafin.com](http://www.biosafin.com) in order to locate the calendar of Courses available.








## Patient Information

For correct Patient information about implant related dental care:

[www.implantologiawinsix.com](http://www.implantologiawinsix.com)



# Summary Table: WINSIX Implants Indications

WINSIX IMPLANTS INDICATIONS DIAGRAM												
Type of Implant	Traditional bone level implant protocol during the osseointegration phase	Immediate post-extraction	Transmucosal Protocol	Platform Switching	Single Tooth Aesthetics	Aesthetics	Single tooth with roots converging apically to adjacent teeth	Direct prosthesis onto implant	Immediate Loading	Just on 4/6	CAB	Multi-unit screw-retained prosthesis
<b>KAPPA line</b>												
<b>K</b>												
	...	...	•	... (not for 3,3)	...	...	-	•	...	-	-	...
<b>KT</b>												
	•	••	...	...	•	•	-	...	...	-	-	...
<b>KE</b>												
	-	...	...	-	...	...	-	...	...	-	-	... (Flat Line)
<b>Torque Type Line</b>												
<b>TTi</b>												
	...	...	•	...	...	...	...	•	...	•	-	...
<b>TTx</b>												
	...	...	•	... (not for 3,3)	-	-	-	...	...	...	...	...
<b>TTc</b>												
	...	...	•	...	...	...	...	-	...	...	...	...
<b>TTcm</b>												
	...	...	...	...	...	...	...	-	...	-	-	-


**KEY**

- Particularly suitable
- Average suitability
- Suitable
- Not suitable

... WINSIX Implant System

# WINSIX Implants: diameters and lengths

The product colour coding system is obtained via a surface anodization process that does not fade on contact with the saliva. It is then possible to quickly identify the diameter of the product accessory when the implant is opened.



## KAPPA line

INTERNAL HEX CONNECTION	Ø 3.3	Ø 3.8	Ø 4.5	Ø 5.2	Ø 5.9
<b>K</b>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	-
	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	-
	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	-
	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	-
<b>KT</b>	<i>h 7</i>	<i>h 7</i>	<i>h 7</i>	<i>h 7</i>	<i>h 7</i>
	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>
	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>
	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>
	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	-
<b>KT MACHINED</b>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	-	-
	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	-	-
	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	-	-
	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	-	-
<b>KE</b>	-	<i>h 6</i>	<i>h 6</i>	<i>h 6</i>	-
	-	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	-
	-	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	-
	-	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	-

## Torque Type Line

INTERNAL HEX CONNECTION	Ø 3.3	Ø 3.8	Ø 4.5	Ø 5.2	Ø 5.9
<b>TORQUE TYPE I</b>	-	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	-
	-	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	-
	-	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	-
	-	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	-
	-	<i>h 18</i>	<i>h 18</i>	-	-

EXTERNAL HEX CONNECTION	Ø 3.3	Ø 3.8	Ø 4.5	Ø 5.2	Ø 5.9
<b>TORQUE TYPE X</b>	-	<i>h 6</i>	<i>h 6</i>	<i>h 6</i>	<i>h 6</i>
	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>
	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>
	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>
	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	-
		<i>h 18</i>	<i>h 18</i>	-	-

CONICAL CONNECTION	Ø 3.3	Ø 3.8	Ø 4.5	Ø 5.2	Ø 5.9
<b>TORQUE TYPE C</b>	-	<i>h 6</i>	<i>h 6</i>	<i>h 6</i>	<i>h 6</i>
	-	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>
	-	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>
	-	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>
	-	<i>h 15</i>	<i>h 15</i>	<i>h 15</i>	-

CONOMETRIC CONNECTION	Ø 3.3	Ø 3.8	Ø 4.5	Ø 5.2	Ø 5.9
<b>TORQUE TYPE CM</b>	-	<i>h 6</i>	<i>h 6</i>	<i>h 6</i>	<i>h 6</i>
	-	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>	<i>h 9</i>
	-	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>	<i>h 11</i>
	-	<i>h 13</i>	<i>h 13</i>	<i>h 13</i>	-


# WINSIX Micro Rough Surface

Over the years, Research and clinical use, the key to the success of WINSIX implants, have proven to be of unique value.


## Excellent osteoconductive properties of the MRS surface

WINSIX Micro Rough Surface presents a rough morphology, obtained by means of the subtractive process of sandblasting and acid-etching directly onto the titanium implant surface. The special three-dimensional geometry of its microstructure allows a reduction in the time required for osseointegration, through a process of osteoconduction promoted by the distinctive scaffold characteristics of this surface.

The regularity of the roughness and the characteristics of MRS promote a greater bone-to-implant contact surface, improving the quality and quantity of bone around the implant and ensuring optimum osseointegration.



MAG 721x: regularity of the MRS surface



Bone showing the active remodelling phase of the MRS surface




The WINSIX Micro Rough Surface, with its special morphology and production Quality - involving 100% screening of the components - is applied to every Implant from the WINSIX Implant System. Through this technology, BIOSAF IN offers a unique value to Dental Surgeons so that they may operate with all the peace of mind they require.

# Implant Collar Types and Heights

The variety of collar type and height offered by WINSIX allows the choice of the most suitable implant for the optimum management of hard and soft tissues.

K




polished collar  
h. 0,3 mm.

KT



polished collar  
h. 1,3 mm.

KE



polished collar  
h. 3,0 mm.

TTi



polished collar  
h. 0,7 mm.

TTx



polished collar  
h. 0,7 mm.

TTc



treated collar

TTcm



treated collar

### Bibliographic References:

• Macroscopic and Microscopic evaluation of a new implant design supporting immediately loaded full arch rehabilitation - S. Tetè, V. Zizzari, A. De Carlo, B. Sinjari, E.F. Gherlone - Annali di Stomatologia Vol.3, no.2 - 2012

• Influence of Novel Nano - Titanium Implant Surface on Human Osteoblast Behaviour and Growth - S. Tetè, F. Mastrangelo, R. Quaresima, R. Vinci, G. Sammartino, L. Stuppia, Enrico F. Gherlone - Implant Dentistry Vol.19, no.26 - 2010

• Isolation of osteogenic progenitors from human amniotic fluid using a single step culture protocol - I. Antonucci, I. Iezzi, E. Morizio, F. Mastrangelo, A. Pantalone, M. Mattioli Belmonte, A. Gigante, V. Salini, G. Calabrese, S. Tetè, G. Palka, L. Stuppia - SILENCE a Journal of RNA Regulation - 2009

• In vitro behaviour onto different titanium surfaces of osteoblast-like cells obtained from human dental pulp - S. Tetè, F. Mastrangelo, V. Zizzari, G. D'Apolito, N. Fiorentino, U. Desiato, M.T. Sberna, R. Quaresima, L. Stuppia, R. Vinci, E.F. Gherlone - Atti del 7th Annual Meeting of ISSCR International Society of Stem Cell Research, Barcelona, July 2009

• Novel Protocol of Osteogenic differentiation from amniotic fluid cells S. Tetè, F. Mastrangelo, M. Tranasi, V. Zizzari, I. Antonucci, G. D'Apolito, T. Marchese, R. Vinci, L. Stuppia, E.F. Gherlone - Atti del 7th Annual Meeting of ISSCR International Society of Stem Cell Research, Barcelona, July 2009

• Interfaccia osso - Implant nei differenti tipi di carico degli impianti dentali S. Tetè, G. D'Apolito, F. Mastrangelo, R. Vinci, E.F. Gherlone - Atti del III Expo di Autunno, Universities of Lombardy, 27-28 November 2009

• Valutazione della capacità osteogenetica di hafscs ottenute da liquido amniotico - S. Tetè, U. Di Tore, V. Zizzari, L. Stuppia, F. Zarone, E.F. Gherlone - Atti del III Expo di Autunno, Universities of Lombardy, 27 - 28 November 2009




# WINSIX Implant Connections Prospectus

The WINSIX System offers 4 different types of connection (internal, external, conical, conometric), in a single dimension in the different implant measurements.

The same instruments and prosthetic components can be used with any type of implant. Ergonomics and system management are optimised and the work of the operators simplified.

# Stable Implant - Abutment fit



## FREE LOCK INTERNAL HEX

- In use since 1995
- It provides great connection stability thanks to the hexagon and screw dimensions.

## EXTERNAL HEX

- 2.7 mm hexagonal wrench. - h. 0.7 mm.
- Particularly suitable for Just on 4/6 Techniques

## INTERNAL CONICAL

- Conicity 3°
- With antirotational hex for the repositioning of prosthetic abutments and for implant screw-retention
- Activated with a low torque value (20 N/cm)
- Low torque value of conicity activation can be very useful in the case of regenerative techniques where bone density is poor
- Platform Switching incorporated

## INTERNAL MORSE TAPER

- Conicity 1.5°
- Percussion activated
- No tightening screw for the prosthetic abutment
- Platform Switching incorporated



# Implant Packaging

The implants are provided in a sterile environment, in a transparent ampoule within a colour coded titanium container. The sterile ampoule is thermo sealed inside a blister pack.

The direct engagement of the implant to its holder is a specific feature of the WINSIX implant. In this way, it:

- has a reduced vertical dimension, facilitating insertion into the oral cavity by means of the insertion instruments.
- never accidentally comes into contact with any materials other than titanium, thus avoiding possible contamination.

- allows the implant to be engaged by the specific mounters, without further manoeuvring during software-assisted implantology.

THE CARDBOARD BOX EXTERIOR reports all the necessary information for the immediate identification of the product (its colour identifies the type of implant it contains), as well as the indicators required by law, in accordance with the standards that regulate medical devices.

The packaging adequately preserves the product, allowing easy storage and immediate visual identification thanks to the colour code that is clearly visible on the box exterior.



the colour of the box identifies the type of implant it contains

K	KT	KTM	KE	TTx	TTi	TTc	TTcm
---	----	-----	----	-----	-----	-----	------

## Inside the Box:

- Information leaflet
- Label with stickers for accurate documental management of Patient data:
  - 1 sticker for the clinical record
  - 1 sticker for communications with the laboratory
  - 1 sticker for the Patient's Personal Record Card



The closure cap supplied with KT implants differs from the implant colour code because it corresponds to the larger prosthesis diameter.



## On the Ampoule:

- Product Code
- Production batch

## On the Blister pack:

- Product Code
- Production batch
- Sterilisation batch
- Sterilisation Date
- Expiry Date



## Colour Code of implant diameters



## Opening the Packaging

The non-sterile assistant opens, first the box, then the blister pack and drops the ampoule onto the surgical tray without touching it.

The sterile operator opens the ampoule by lifting the cap to which the implant closure cap is attached.

The implant is housed inside the sterile ampoule. To remove it, use the appropriate tools without turning it over.

## Extraction of the implant

The implant is directly disengaged by the operator from the centre of the titanium holder using the handpiece, manually or by means of a ratchet, without interrupting the sterile chain.

# Implants


## KAPPA Line

# KAPPA line

3 DIFFERENT IMPLANT COLLARS  
THE SAME IMPLANT BODY AND SAME CONNECTION  
Internal HEX



**Varying DEPTHS**  
The THREADS alter gradually in depth to promote progressive horizontal micro-expansion.



## macromorphology

### Variable GEOMETRY

The THREADS gradually alter in form from quadrilateral to triangular to promote vertical micro-expansion.

### Micro Rough Surface

With rough morphology, obtained by means of the subtractive process applied to the implant surface.

### DISCHARGE GROOVES

Wide and deep for the deposit of bone fragments and blood clot formation during implant insertion.



### OPTIMUM STABILITY

in any kind of bone, thanks to threads of variable geometry and varying depths. Also, ideal for immediate loading.

K

polished collar  
h. 0,3 mm.




### Crestal module

With microgrooves for greater bone stability in the coronal area and subsequently, excellent aesthetic results.

KT

polished collar  
h. 1,3 mm.




### Crestal module

with platform that switches to the consecutive prosthetic diameter. Ideal for direct application of the prosthesis onto the implant with multi-unit screw-retained prosthesis. It combines optimum functionality and cost reduction. Useful in the distal area to reproduce a prosthetic crown anatomy similar to the natural tooth anatomy of molars.

KE

polished collar  
h. 3,0 mm.




### Crestal module ideal for:



- prosthetic application with direct fastening to the implant collar which enhances the aesthetic aspects to the highest level.
- direct application of the prosthesis onto the implant without prosthetic transmucosal abutment.

It combines optimum functionality and cost reduction with the utmost regard for the biological width of soft tissues.

# ...K Implant



- Placement:** bone level
- Shape:** cylindrical with truncated-cone body
- Collar:** aesthetic h. 0.3 mm with microgrooves for greater bone stability in the coronal area and subsequent excellent aesthetic outcome.
- Connection:** Free Lock internal hex
- Implant body:** self-tapping thread, of varying depth and thickness, and variable geometry to modulate primary stability in surgical phases. Optimum stability in any type of bone, excellent for immediate loading.
- Surface:** WINSIX Micro Rough Surface (MRS)
- Apex:** semi-spherical, with wide and deep discharge grooves (1/3 of the implant body) for bone fragment collection ideal for maxillary sinus lift surgery.



internal

Product Code

Product

**K IMPLANTS**  
of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h. 0.3 mm, closure cap included (code VT..)



33009K/MRS  
33011K/MRS  
33013K/MRS  
33015K/MRS

- 3.3 mm. x h. 9 mm.
- 3.3 mm. x h.11 mm.
- 3.3 mm. x h.13 mm.
- 3.3 mm. x h.15 mm.

38009K/MRS  
38011K/MRS  
38013K/MRS  
38015K/MRS

- 3.8 mm. x h. 9 mm.
- 3.8 mm. x h.11 mm.
- 3.8 mm. x h.13 mm.
- 3.8 mm. x h.15 mm.

45009K/MRS  
45011K/MRS  
45013K/MRS  
45015K/MRS


- 4.5 mm. x h. 9 mm.
- 4.5 mm. x h.11 mm.
- 4.5 mm. x h.13 mm.
- 4.5 mm. x h.15 mm.

52009K/MRS  
52011K/MRS  
52013K/MRS  
52015K/MRS


- 5.2 mm. x h. 9 mm.
- 5.2 mm. x h.11 mm.
- 5.2 mm. x h.13 mm.
- 5.2 mm. x h.15 mm.



# ...KT Implant



- Placement:** transmucosal
- Shape:** cylindrical with truncated-cone body
- Collar:** flared, polished, transmucosal h. 1.3 mm. with prosthetic platform with a diameter wider than the endosseous implant diameter.
- ADVANTAGES:**
- Optional benefit of use as bone level implant.
  - For better prosthetic management in case of reduced bone depth.
- Connection:** Free Lock internal hex
- Implant body:** self-tapping thread, of varying depth and thickness, and variable geometry to modulate primary stability in surgical phases. Optimum stability in any type of bone, excellent for immediate loading.
- Surface:** WINSIX Micro Rough Surface (MRS)
- Apex:** semi-spherical, with wide and deep discharge grooves (1/3 of the implant body) for bone fragment collection ideal for maxillary sinus lift surgery.



GH Cylindrical closure cap.2 mm.

Alternatively, the KT implant can be submerged to bone level using the flat closure cap which is normally fitted with other implants and can be purchased separately.



## Product Code

## Product

33007KT/MRS  
33009KT/MRS  
33011KT/MRS  
33013KT/MRS  
33015KT/MRS

**KT IMPLANTS**  
of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h.1.3 mm., healing caps included (code. VG3802)

- 3.3 mm. x h. 7 mm.
- 3.3 mm. x h. 9 mm.
- 3.3 mm. x h. 11 mm.
- 3.3 mm. x h. 13 mm.
- 3.3 mm. x h. 15 mm.

38007KT/MRS  
38009KT/MRS  
38011KT/MRS  
38013KT/MRS  
38015KT/MRS



**KT IMPLANTS**  
of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h.1.3 mm., healing caps included (code VG4502)

- 3.8 mm. x h. 7 mm.
- 3.8 mm. x h. 9 mm.
- 3.8 mm. x h. 11 mm.
- 3.8 mm. x h. 13 mm.
- 3.8 mm. x h. 15 mm.



# ...KT Implant

# ...KT Machined Implant



Product Code

Product

Product Code

Product

### KT IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h.1.3 mm., healing caps included (code VG5202)

- 4.5 mm. x h. 7 mm.
- 4.5 mm. x h. 9 mm.
- 4.5 mm. x h.11 mm.
- 4.5 mm. x h.13 mm.
- 4.5 mm. x h.15 mm.

### KT IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h.1.3 mm., healing caps included (code VG5902)

- 5.2 mm. x h. 7 mm.
- 5.2 mm. x h. 9 mm.
- 5.2 mm. x h.11 mm.
- 5.2 mm. x h.13 mm.
- 5.2 mm. x h.15 mm.

### KT IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h.1.3 mm., healing caps included (code VG6502)

- 5.9 mm. x h. 7 mm.
- 5.9 mm. x h. 9 mm.
- 5.9 mm. x h.11 mm.
- 5.9 mm. x h.13 mm.

### KTM IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, Machined surface, closure cap included (code VT38)

- 3.3 mm. x h. 9 mm.
- 3.3 mm. x h. 11 mm.
- 3.3 mm. x h. 13 mm.
- 3.3 mm. x h. 15 mm.

### KTM IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, Machined surface, closure cap included (code VT45)

- 3.8 mm. x h. 9 mm.
- 3.8 mm. x h. 11 mm.
- 3.8 mm. x h. 13 mm.
- 3.8 mm. x h. 15 mm.

### KTM IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, Machined surface, closure cap included (code VT52)

- 4.5 mm. x h. 9 mm.
- 4.5 mm. x h.11 mm.
- 4.5 mm. x h.13 mm.
- 4.5 mm. x h.15 mm.



45007KT/MRS  
45009KT/MRS  
45011KT/MRS  
45013KT/MRS  
45015KT/MRS



52007KT/MRS  
52009KT/MRS  
52011KT/MRS  
52013KT/MRS  
52015KT/MRS



59007KT/MRS  
59009KT/MRS  
59011KT/MRS  
59013KT/MRS



33009KTM  
33011KTM  
33013KTM  
33015KTM




38009KTM  
38011KTM  
38013KTM  
38015KTM



45009KTM  
45011KTM  
45013KTM  
45015KTM

# KE Implant

 internal



- Placement:** transmucosal
- Shape:** cylindrical with truncated-cone body
- Collar:** aesthetic h. 3.0 mm., allows placement of the prosthetic crown directly onto the collar. Excellent aesthetic result, thanks to the absence of a micro gap between the abutment and the implant.
- Connection:** Free Lock internal hex
- Implant body:** self-tapping thread, of varying depth and thickness, and variable geometry to modulate primary stability in surgical phases. Optimum stability in any type of bone, excellent for immediate loading.
- Surface:** WINSIX Micro Rough Surface (MRS)
- Apex:** semi-spherical, with wide and deep discharge grooves (1/3 of the implant body) for the collection of bone fragments.



## Product Code

## Product

### KE IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h. 3,0 mm., closure cap included (cod VT33)

38006KE/MRS  
38009KE/MRS  
38011KE/MRS  
38013KE/MRS

- 3.8 mm. x h. 6 mm.
- 3.8 mm. x h. 9 mm.
- 3.8 mm. x h.11 mm.
- 3.8 mm. x h.13 mm.

### KE IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h. 3,0 mm., closure cap included (code VT38)

45006KE/MRS  
45009KE/MRS  
45011KE/MRS  
45013KE/MRS

- 4.5 mm. x h. 6 mm.
- 4.5 mm. x h. 9 mm.
- 4.5 mm. x h.11 mm.
- 4.5 mm. x h.13 mm.

### KE IMPLANTS

of cylindrical shape with truncated-cone body self-tapping, MRS surface, polished collar h. 3,0 mm., closure cap included (code VT45)

52006KE/MRS  
52009KE/MRS  
52011KE/MRS  
52013KE/MRS


- 5.2 mm. x h. 6 mm.
- 5.2 mm. x h. 9 mm.
- 5.2 mm. x h.11 mm.
- 5.2 mm. x h.13 mm.

# KAPPA Line


## Soft tissue management

The KAPPA Line implants share the same implant body and are differentiated by type of collar, offering optimum soft tissue management.


### Implant K



### Implant KT



### Implant KE



## Colour Code

IMPLANT K		
NOMINAL	ENDOSSEOUS Ø	PROSTHETIC PLATFORM Ø
Ø 3.3	● 3.5	● 3.3
Ø 3.8	● 4.0	● 3.8
Ø 4.5	● 4.7	● 4.5
Ø 5.2	● 5.4	● 5.2

**SPECIAL FEATURES OF THE K IMPLANT:**  
 With microgrooves for greater bone stability in the coronal area and subsequently, excellent aesthetic results.  
 • Semi-spherical apex particularly suitable for maxillary sinus lift, to minimise the possibility of sinus membrane laceration.

IMPLANT KT		
NOMINAL	ENDOSSEOUS Ø	PROSTHETIC PLATFORM Ø
Ø 3.3	● 3.5	● 3.8
Ø 3.8	● 4.0	● 4.5
Ø 4.5	● 4.7	● 5.2
Ø 5.2	● 5.4	● 5.9
Ø 5.9	● 6.1	● 6.5

**SPECIAL FEATURES OF THE KT IMPLANT:**  
 • provides the possibility of customised management of peri-implant tissue during the healing phase. In fact, the Closure cap, which is cylindrical in form with h. 2 mm. carries out a transgingival function.  
 • possibility of using a wider prosthetic Ø even in cases where there is less bone availability.  
 • optimum harmony between the biological width of the soft tissue and the aesthetic result.  
 • platform switching possible from endosseous Ø of 3.3 mm.

IMPLANT KE		
NOMINAL	ENDOSSEOUS Ø	PROSTHETIC PLATFORM Ø
Ø 3.8	● 4.0	● 3.3
Ø 4.5	● 4.7	● 3.8
Ø 5.2	● 5.4	● 4.5

**SPECIAL FEATURES OF THE K IMPLANT:**  
 • Possible application of the prosthetic crown directly onto the implant collar, eliminating any crevices at the crestal level and between the abutment and the implant.  
 • Possible application of the multi-unit screw-retained prostheses, directly onto the implant, perfectly combining the precepts of biological width, soft tissues and cost reduction.

# ...Implants

TORQUE TYPE LINE



# TORQUE TYPE Line

**TTi**  
internal hex


**just one implant body**

specific macromorphology to achieve maximum implant stability, also ideal for immediate loading.

**TTx**  
external hex

**POLISHED COLLAR**

polished collar truncated-cone shape for optimum soft tissue management



**External Hex**  
○ constant 2.7 mm height 0.7 mm

**BIO-ENGAGING THREADS**

**Double-threaded and double-step threads** for easy implant placement with half the turns required. The groove at the bottom of the thread decompresses the bone by dissipating the forces and facilitating blood clot deposit. At the same time, it increases the implant surface by facilitating cell neoformation.

**Multifunctional tapered apex**  
Ø undersized from 1.3 to 1.8 mm compared to the Ø of the implant. Heavily tapered in order to obtain an osteotomic effect and facilitate tilted implant placement even in cases of reduced bone availability, according to Teeth Just on 4/6 Techniques. It modulates primary stability in relation to the bone type. Ideal for:  
- apically convergent roots  
- Apically convergent implants



**TAPERED APEX**

**120° Discharge grooves**  
With osteotomic bone expansion and decompression function. Wide and deep for the harvesting of bone fragments and deposit of blood clots during implant insertion.

**IDEAL FOR JUST ON 4/6 TECHNIQUES**



**OPTIMUM STABILITY**

in any type of bone, thanks to undersized apex of 1.8 mm. and the osteotomic threads. Also, ideal for immediate loading.

**TTi**




**Internal hex connection**  
Free Lock Connection (equal to KAPPA Line).



**TTx**



**External hex connection**  
Ideal for Just on 4/6 technique. Designed for restorations with CAB device



**TTc**




**Internal conical connection with antirotational hex** (conicity 3°)  
Suitable for prosthetic superstructures with conical activation.



**TTcm**



**Conometric Connection** (conicity 1.5°)  
Suitable for aesthetic solutions. No transverse screws.





**Placement:** bone level

**Shape:** cylindrical with tapered apex


**Collar:** polished h. 0.7 mm. truncated-cone shape

**Connection:** Free Lock internal hex

**Implant body:** double-threaded and double-step threads for facilitated implant placement with half the turns required.

**Surface:** WINSIX Micro Rough Surface (MRS)

**Apex:** multifunctional tapered, undersized  $\varnothing$  of 1.8 mm compared to the  $\varnothing$  of the implant.  
Ideal for:  
- apically convergent roots  
- apically convergent implants



## Product Code

## Product

### TTi IMPLANTS

cylindrical shape with tapered apex  
Free Lock internal connection  
self-tapping, MRS surface polished collar h. 0.7 mm,  
closure cap included (cod VT ..)



38009TTi/MRS  
38011TTi/MRS  
38013TTi/MRS  
38015TTi/MRS  
38018TTi/MRS

● 3.8 mm. x h. 9 mm.  
● 3.8 mm. x h. 11 mm.  
● 3.8 mm. x h. 13 mm.  
● 3.8 mm. x h. 15 mm.  
● 3.8 mm. x h. 18 mm.



45009TTi/MRS  
45011TTi/MRS  
45013TTi/MRS  
45015TTi/MRS  
45018TTi/MRS

● 4.5 mm. x h. 9 mm.  
● 4.5 mm. x h. 11 mm.  
● 4.5 mm. x h. 13 mm.  
● 4.5 mm. x h. 15 mm.  
● 4.5 mm. x h. 18 mm.



52009TTi/MRS  
52011TTi/MRS  
52013TTi/MRS  
52015TTi/MRS

● 5.2 mm. x h. 9 mm.  
● 5.2 mm. x h. 11 mm.  
● 5.2 mm. x h. 13 mm.  
● 5.2 mm. x h. 15 mm.


# TTx Implant





- Placement:** bone level
- Shape:** cylindrical with tapered apex
- Collar:** polished h. 0.7 mm. truncated-cone shape
- Connection:** external hex
- Implant body:** double-threaded and double-step threads for facilitated implant placement with half the turns required.
- Surface:** WINSIX Micro Rough Surface (MRS)
- Apex:** multifunctional tapered, undersized  $\varnothing$  from 1.3 to 1.8 mm compared to the  $\varnothing$  of the implant.  
Ideal for:
  - apically convergent roots
  - apically convergent implants



External connection of the TORQUE TYPE x implant with a constant  $\varnothing$  of 2.7 mm, height 0.7 mm



TORQUE TYPE X IMPLANTS with external hex connection, and constant hex dimension for the prosthetic connection.



## Product Code

## Product

### TTx IMPLANTS

of cylindrical shape with tapered apex, self-tapping external hex connection, MRS surface, polished collar h. 0.7 mm, closure cap included (code VT..x)



33009TTx/MRS  
33011TTx/MRS  
33013TTx/MRS  
33015TTx/MRS

- 3.3 mm. x h. 9 mm.
- 3.3 mm. x h. 11 mm.
- 3.3 mm. x h. 13 mm.
- 3.3 mm. x h. 15 mm.



38006TTx/MRS  
38009TTx/MRS  
38011TTx/MRS  
38013TTx/MRS  
38015TTx/MRS  
38018TTx/MRS

- 3.8 mm. x h. 6 mm.
- 3.8 mm. x h. 9 mm.
- 3.8 mm. x h. 11 mm.
- 3.8 mm. x h. 13 mm.
- 3.8 mm. x h. 15 mm.
- 3.8 mm. x h. 18 mm.



45006TTx/MRS  
45009TTx/MRS  
45011TTx/MRS  
45013TTx/MRS  
45015TTx/MRS  
45018TTx/MRS

- 4.5 mm. x h. 6 mm.
- 4.5 mm. x h. 9 mm.
- 4.5 mm. x h. 11 mm.
- 4.5 mm. x h. 13 mm.
- 4.5 mm. x h. 15 mm.
- 4.5 mm. x h. 18 mm.



52006TTx/MRS  
52009TTx/MRS  
52011TTx/MRS  
52013TTx/MRS  
52015TTx/MRS

- 5.2 mm. x h. 6 mm.
- 5.2 mm. x h. 9 mm.
- 5.2 mm. x h. 11 mm.
- 5.2 mm. x h. 13 mm.
- 5.2 mm. x h. 15 mm.



59006TTx/MRS  
59009TTx/MRS  
59011TTx/MRS  
59013TTx/MRS

- 5.9 mm. x h. 6 mm.
- 5.9 mm. x h. 9 mm.
- 5.9 mm. x h. 11 mm.
- 5.9 mm. x h. 13 mm.

TTx height 6 mm.



### TORQUE TYPE SHORT IMPLANT

- used in the event of minimally invasive interventions
- under the maxillary sinuses
- above the mandibular nerve

# TTc Implant

 conical



- Placement:** bone level
- Shape:** cylindrical with tapered apex
- Collar:** fully treated, truncated-cone shape, MRS surface, optimum bone stability in coronal area, excellent aesthetic result
- Connection:** internal cone with antirotational hex (conicity 3°) and incorporated platform switching
- Implant body:** double-threaded and double-step threads for easy implant placement with half the turns required.
- Surface:** WINSIX Micro Rough Surface (MRS)
- Apex:** multifunctional tapered, undersized  $\varnothing$  of 1.8 mm compared to the  $\varnothing$  of the implant.  
Ideal for:
  - apically convergent roots
  - apically convergent implants



## Product Code

## Product


### TTc IMPLANTS

Cylindrical shape with tapered apex, conical connection, self-tapping, MRS surface, fully treated collar including closure cap (code VT..c)

38006TTc/MRS	● 3.8 mm. x h. 6 mm.
38009TTc/MRS	● 3.8 mm. x h. 9 mm.
38011TTc/MRS	● 3.8 mm. x h.11 mm.
38013TTc/MRS	● 3.8 mm. x h.13 mm.
38015TTc/MRS	● 3.8 mm. x h.15 mm.
45006TTc/MRS	● 4.5 mm. x h. 6 mm.
45009TTc/MRS	● 4.5 mm. x h. 9 mm.
45011TTc/MRS	● 4.5 mm. x h.11 mm.
45013TTc/MRS	● 4.5 mm. x h.13 mm.
45015TTc/MRS	● 4.5 mm. x h.15 mm.
52006TTc/MRS	● 5.2 mm. x h. 6 mm.
52009TTc/MRS	● 5.2 mm. x h. 9 mm.
52011TTc/MRS	● 5.2 mm. x h.11 mm.
52013TTc/MRS	● 5.2 mm. x h.13 mm.
52015TTc/MRS	● 5.2 mm. x h.15 mm.
59006TTc/MRS	● 5.9 mm. x h. 6 mm.
59009TTc/MRS	● 5.9 mm. x h. 9 mm.
59011TTc/MRS	● 5.9 mm. x h.11 mm.
59013TTc/MRS	● 5.9 mm. x h.13 mm.



# TTcm Implant



- Placement:** bone level
- Shape:** cylindrical with tapered apex
- Collar:** truncated-cone shape, fully treated, MRS surface, optimum bone stability in coronal area, excellent aesthetic result
- Connection:** Morse taper, conicity 1.5°
- Implant body:** double-threaded and double-step threads for facilitated implant placement with half the turns required.
- Surface:** WINSIX Micro Rough Surface (MRS)
- Apex:** multifunctional tapered, undersized  $\varnothing$  of 1.8 mm compared to the  $\varnothing$  of the implant.  
Ideal for:
  - apically convergent roots
  - apically convergent implants



## Product Code

## Product

### TTcm IMPLANTS

Cylindrical shape with tapered apex, morse taper connection, self-tapping, MRS surface, fully treated collar including closure cap (code VT..cm)

38006TTcm/MRS	● 3.8 mm. x h. 6 mm.
38009TTcm/MRS	● 3.8 mm. x h. 9 mm.
38011TTcm/MRS	● 3.8 mm. x h. 11 mm.
38013TTcm/MRS	● 3.8 mm. x h. 13 mm.
45006TTcm/MRS	● 4.5 mm. x h. 6 mm.
45009TTcm/MRS	● 4.5 mm. x h. 9 mm.
45011TTcm/MRS	● 4.5 mm. x h. 11 mm.
45013TTcm/MRS	● 4.5 mm. x h. 13 mm.
52006TTcm/MRS	● 5.2 mm. x h. 6 mm.
52009TTcm/MRS	● 5.2 mm. x h. 9 mm.
52011TTcm/MRS	● 5.2 mm. x h. 11 mm.
52013TTcm/MRS	● 5.2 mm. x h. 13 mm.
59006TTcm/MRS	● 5.9 mm. x h. 6 mm.
59009TTcm/MRS	● 5.9 mm. x h. 9 mm.
59011TTcm/MRS	● 5.9 mm. x h. 11 mm.



# ...Surgery

INSTRUMENTS

ST16 SURGICAL TRAY




- STP6N drill stop for 6 mm implants TT1S-TT2S-TT3S
- STP7N drill stop for 7 mm implants TT1S-TT2S-TT3S
- STP9N drill stop for 9 mm implants TT1S-TT2S-TT3S
- STP11N drill stop for 11 mm implants TT1S-TT2S-TT3S
  
- STP6R drill stop for 6 mm implants TT4S-TT5S-TT6S
- STP7R drill stop for 7 mm implants TT4S-TT5S-TT6S
- STP9R drill stop for 9 mm implants TT4S-TT5S-TT6S
- STP11R drill stop for 11 mm implants TT4S-TT5S-TT6S
  
- STP6W drill stop for 6 mm implants TT7S-TT8S-TT9-TT10
- STP7W drill stop for 7 mm implants TT7S-TT8S-TT9-TT10
- STP9W drill stop for 9 mm implants TT7S-TT8S-TT9-TT10
- STP11W drill stop for 11 mm implants TT7S-TT8S-TT9-TT10

60000 Ratchet

KMP adapter for ratchet screwdriver

65000K dynamometric ratchet 10 - 35 N/cm





## Available materials

Product Code	Description	N.
TTK	Precision drill	1
TT1S	Short drill Ø: 2.0 mm - L: 6-7-9-11-13 mm	1
TT2S	Short drill Ø: 2.2/2.6 mm - L: 6-7-9-11-13 mm	1
TT3S	Short drill Ø: 2.6/3.0 mm - L: 6-7-9-11-13 mm	1
TT4S	Short drill Ø: 3.0/3.4 mm - L: 6-7-9-11-13 mm	1
TT5S	Short drill Ø: 3.4/3.8 mm - L: 6-7-9-11-13 mm	1
TT6S	Short drill Ø: 3.8/4.2 mm - L: 6-7-9-11-13 mm	1
TT7S	Short drill Ø: 4.2/4.6 mm - L: 6-7-9-11-13 mm	1
TT8S	Short drill Ø: 4.6/5.0 mm - L: 6-7-9-11-13 mm	1
TT9	Short drill Ø: 5.0/5.4 mm - L: 6-7-9-11-13 mm	1
TT10	Short drill Ø: 5.4/5.8 mm - L: 6-7-9-11-13 mm	1
TT1L	Long drill Ø: 2.0 mm - L: 6-7-9-11-13-15-18 mm	1
TT2L	Long drill Ø: 2.2/2.6 mm - L: 6-7-9-11-13-15-18 mm	1
TT3L	Long drill Ø: 2.6/3.0 mm - L: 6-7-9-11-13-15-18 mm	1
TT4L	Long drill Ø: 3.0/3.4 mm - L: 6-7-9-11-13-15-18 mm	1
TT5L	Long drill Ø: 3.4/3.8 mm - L: 6-7-9-11-13-15-18 mm	1
TT6L	Long drill Ø: 3.8/4.2 mm - L: 6-7-9-11-13-15-18 mm	1
TT7L	Long drill Ø: 4.2/4.6 mm - L: 6-7-9-11-13-15-18 mm	1
TT8L	Long drill Ø: 4.6/5.0 mm - L: 6-7-9-11-13-15-18 mm	1
99000	Drill Extension Shaft	1
BM33	Bone mill for EA Ø 3.3 mm	1
BM38	Bone mill for EA Ø 3.8 mm	1
BM45	Bone mill for EA Ø 4.5 mm	1
BM52	Bone mill for EA Ø 5.2 mm	1
BM59	Bone mill for EA Ø 5.9 mm	1
CD33	Cortical Drill Ø 3.3 mm	1
CD38	Cortical Drill Ø 3.8 mm	1
CD45	Cortical Drill Ø 4.5 mm	1
CD52	Cortical Drill Ø 5.2 mm	1
CD59	Cortical Drill Ø 5.9 mm	1
DI	Direction Indicator	4
DPI	Depth Indicator	1
62SI	Short Screwdriver for internal hex connection Implants	1
62MI	Medium Screwdriver for internal hex connection Implants	1
62LI	Long Screwdriver for internal hex connection Implants	1
99SI	Short Handpiece Screwdriver for internal hex connection Implants	1
99MI	Medium Handpiece Screwdriver for internal hex connection Implants	1
99LI	Long Handpiece Screwdriver for internal hex connection Implants	1
62SX	Short Screwdriver for external hex connection Implants	1
62MX	Medium Screwdriver for external hex connection Implants	1
62LX	Long Screwdriver for external hex connection Implants	1











## for the Surgical Kit









surgery


Product Code	Description	N.
99SX	Short Handpiece Screwdriver for external hex connection Implants	1
99MX	Medium Handpiece Screwdriver for external hex connection Implants	1
99LX	Long Handpiece Screwdriver for external hex connection Implants	1
62Sc	Short Screwdriver for conical connection Implants	1
62Mc	Medium Screwdriver for conical connection Implants	1
62Lc	Long Screwdriver for conical connection Implants	1
99Sc	Short Handpiece Screwdriver for conical connection Implants	1
99Mc	Medium Handpiece Screwdriver for conical connection Implants	1
99Lc	Long Handpiece Screwdriver for conical connection Implants	1
62Scm	Short Screwdriver for conometric connection Implants	1
62Mcm	Medium Screwdriver for conometric connection Implants	1
62Lcm	Long screwdriver for conometric connection Implants	1
99Scm	Short Handpiece Screwdriver for morse taper connection Implants	1
99Mcm	Medium Handpiece Screwdriver for morse taper connection Implants	1
99Lcm	Long Handpiece Screwdriver for morse taper connection Implants	1
SDL	Long Manual Screwdriver	1
SDS	Short Manual Screwdriver	1
30000	Prosthetic wrench with handpiece connection	1
SDTL	Long Torx Manual Screwdriver	1
SDTS	Short Torx Manual Screwdriver	1
30000T	Torx Wrench with handpiece connection	1
60000	Ratchet	1
KMP	Adapter for ratchet screwdriver	1
STP6N	6 mm implant drill stops for TT1S-TT2S-TT3S drills	1
STP7N	7 mm implant drill stops for TT1S-TT2S-TT3S drills	1
STP9N	9 mm implant drill stops for TT1S-TT2S-TT3S drills	1
STP11N	11 mm implant drill stops for TT1S-TT2S-TT3S drills	1
STP6R	6 mm implant drill stops for TT4S-TT5S-TT6S drills	1
STP7R	7 mm implant drill stops for TT4S-TT5S-TT6S drills	1
STP9R	9 mm implant drill stops for TT4S-TT5S-TT6S drills	1
STP11R	11 mm implant drill stops for TT4S-TT5S-TT6S drills	1
STP6W	6 mm implant drill stops for TT7S-TT8S-TT9-TT10 drills	1
STP7W	7 mm implant drill stops for TT7S-TT8S-TT9-TT10 drills	1
STP9W	9 mm implant drill stops for TT7S-TT8S-TT9-TT10 drills	1
STP11W	11 mm implant drill stops for TT7S-TT8S-TT9-TT10 drills	1
62Ttc	Manual Extractor for conical connection	1
62EA	Manual screwdriver for straight EA	1
99EA	Straight EA screwdriver for handpiece	1
ST16	Single Surgical Tray	1



## Summary table for the use of Drills KAPPA line

RANGE INDICATING THE CORRECT DRILL TO USE IN THE EVENT OF:		SOFT BONE		MEDIUM BONE			HARD BONE			
		● ● ● ● ●	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●	● ● ● ● ●
<b>Short Drill</b>										
PRECISION DRILL TTK	TT1S	TT2S	TT3S	TT4S	TT5S	TT6S	TT7S	TT8S	TT9	TT10
										
	2.0 mm.	2.2/2.6 mm.	2.6/3.0 mm.	3.0/3.4 mm.	3.4/3.8 mm.	3.8/4.2 mm.	4.2/4.6 mm.	4.6/5.0 mm.	5.0/5.4 mm.	5.4/5.8 mm.

<b>Long Drill</b>										
PRECISION DRILL TTK	TT1L	TT2L	TT3L	TT4L	TT5L	TT6L	TT7L	TT8L		
										
	2.0 mm.	2.2/2.6 mm.	2.6/3.0 mm.	3.0/3.4 mm.	3.4/3.8 mm.	3.8/4.2 mm.	4.2/4.6 mm.	4.6/5.0 mm.		

<b>Cortical Drill</b>					
 cortical drill CD 3.3 mm.	 cortical drill CD 3.8 mm.	 cortical drill CD 4.5 mm.	 cortical drill CD 5.2 mm.	 cortical drill CD 5.9 mm.	Recommended for compact cortical tissue and for all implants
					→

## Summary table for the use of Drills KAPPA line

surgery

SOFT BONE				
Ø Nominal diameter of the implant	Ø Endosseous diameter	Ø Drill diameter	Last Drill Code	Cortical Drill
● 3.3	3.5	2.2/2.6	TT2	NO
● 3.8	4.0	2.6/3.0	TT3	NO
● 4.5	4.7	3.0/3.4	TT4	NO
● 5.2	5.4	3.4/3.8	TT5	NO
● 5.9	6.1	3.8/4.2(4.2/4.6)	TT6(TT7)	NO
MEDIUM BONE				
Ø Nominal diameter of the implant	Ø Endosseous diameter	Ø Drill diameter	Last Drill Code	Cortical Drill
● 3.3	3.5	2.6/3.0	TT3	NO
● 3.8	4.0	3.0/3.4	TT4	NO
● 4.5	4.7	3.4/3.8	TT5	NO
● 5.2	5.4	3.8/4.2	TT6	NO
● 5.9	6.1	4.2/4.6(4.6/5.0)	TT7(TT8)	NO
MEDIUM CORTICAL BONE				
Ø Nominal diameter of the implant	Ø Endosseous diameter	Ø Drill diameter	Last Drill Code	Cortical Drill
● 3.3	3.5	2.6/3.0	TT3	CD33
● 3.8	4.0	3.0/3.4	TT4	CD38
● 4.5	4.7	3.8/4.2	TT6	CD45
● 5.2	5.4	4.6/5.0	TT8	CD52
● 5.9	6.1	5.0/5.4	TT9	CD59
HARD BONE				
Ø Nominal diameter of the implant	Ø Endosseous diameter	Ø Drill diameter	Last Drill Code	Cortical Drill
● 3.3	3.5	3.0/3.4	TT4	
● 3.8	4.0	3.4/3.8	TT5	
● 4.5	4.7	3.8/4.2	TT6	CD45
● 5.2	5.4	4.6/5.0	TT8	
● 5.9	6.1	5.4/5.8	TT10	

Drills TT		Cortical Drills	
Code	Ø	Code	Ø
TT1	2.0		
TT2	2.2/2.6		
TT3	2.6/3.0	CD33	3.0/3.3
TT4	3.0/3.4	CD38	3.4/3.8
TT5	3.4/3.8		
TT6	3.8/4.2	CD45	4.2/4.5
TT7	4.2/4.6		
TT8	4.6/5.0	CD52	5.0/5.2
TT9	5.0/5.4	CD59	5.4/5.9
TT10	5.4/5.8		



## Precision

## Long

Product Code

Product

Product Code

Product




TTK

**PRECISION DRILL**  
in surgical steel, sterilisable,  
cutting edge length 5 mm.

- TT1L
- TT2L
- TT3L
- TT4L
- TT5L
- TT6L
- TT7L
- TT8L


**LONG DRILL**  
in surgical steel, sterilisable,  
of variable length

- Ø 2.0 mm. L. 6-18 mm.
- Ø 2.2/2.6 mm. L. 6-18 mm.
- Ø 2.6/3.0 mm. L. 6-18 mm.
- Ø 3.0/3.4 mm. L. 6-18 mm.
- Ø 3.4/3.8 mm. L. 6-18 mm.
- Ø 3.8/4.2 mm. L. 6-18 mm.
- Ø 4.2/4.6 mm. L. 6-18 mm.
- Ø 4.6/5.0 mm. L. 6-18 mm.



The drill lengths are overestimated by 0.5 mm. compared to the endosseous length of the implants.

## Short



TT Short Drill Depth markings


The drill lengths are overestimated by 0.5 mm. compared to the endosseous length of the implants.

- TT1S
- TT2S
- TT3S
- TT4S
- TT5S
- TT6S
- TT7S
- TT8S
- TT9
- TT10

**SHORT DRILLS**  
in surgical steel, sterilisable, of variable length  
Depth stop of 13 mm. included with the drill

- Ø 2.0 mm. L.6-13 mm.
- Ø 2.2/2.6 mm. L.6-13 mm.
- Ø 2.6/3.0 mm. L.6-13 mm.
- Ø 3.0/3.4 mm. L.6-13 mm.
- Ø 3.4/3.8 mm. L.6-13 mm.
- Ø 3.8/4.2 mm. L.6-13 mm.
- Ø 4.2/4.6 mm. L.6-13 mm.
- Ø 4.6/5.0 mm. L.6-13 mm.
- Ø 5.0/5.4 mm. L.6-13 mm.
- Ø 5.4/5.8 mm. L.6-13 mm.

## Cortical



- CD33
- CD38
- CD45
- CD52
- CD59

**CORTICAL DRILLS**  
in surgical steel, sterilisable, L. 5 mm.

- 3.3 mm. L. 5 mm.
- 3.3 mm. L. 5 mm.
- 4.5 mm. L. 5 mm.
- 5.2 mm. L. 5 mm.
- 5.9 mm. L. 5 mm.

## Product Code

## Product



BM33  
BM38  
BM45  
BM52  
BM59

### BONE MILLS FOR EA

in surgical steel, sterilisable

- **3.3 mm.** associated with code VG3305X - VG3305
- **3.8 mm.** associated with code VG3805X - VG3805
- **4.5 mm.** associated with code VG4505X - VG4505
- **5.2 mm.** associated with code VG5205X - VG5205
- **5.9 mm.** associated with code VG5905X - VG5905

The bone mill must be used together with the straight 5 mm Healing Caps. These are inserted into the implant.

## Product Code

## Product

### IMPLANT DRILL STOPS

TT1S - TT2S - TT3S

in surgical steel, sterilisable

da 6 mm.  
da 7 mm.  
da 9 mm.  
da 11 mm.

### STOPS FOR SHORT DRILLS

TT4S - TT5S - TT6S

in surgical steel, sterilisable


da 6 mm.  
da 7 mm.  
da 9 mm.  
da 11 mm.


### STOPS FOR SHORT DRILLS

TT7S - TT8S - TT9 - TT10

in surgical steel, sterilisable

da 6 mm.  
da 7 mm.  
da 9 mm.  
da 11 mm.





Guided tilting for drills

Product Code

Product

WG46



**WINSIX Just Guide for Just on 4/6**  
guides the inclination of the drill from 0° to 45° by means of a specific laser marking. Can be completely dismantled for washing and sterilisation. To reassemble it use the Short/Long Torx Manual Screwdriver (SDTS-SDTL)

DI



**Direction Indicator**  
usable on both sides, to verify the axis of the surgical site after use of drills with Ø 2.0. Equipped with a central through hole, total length 25 mm.

DPI



**Depth Indicator**  
to verify the preparation depth of the surgical site from 6 to 18 mm., after use of drills with Ø 2.0.

99EA



**Straight EA Screwdriver for handpiece**

62EA



**Ratchet screwdriver for straight EA**

Straight EAs should be tightened to C 30 N/cm

62TTC



**Manual Extractor for conical connection**

Product Code

Product

**VISUAL CARD**

**The Visual Card accurately reproduces** the implant series in their various sizes in the following scales: 1: 1 - 1.1: 1 - 1.2: 2. During the planning phase, this is overlapped by the patient's radiographs which allows assessment of the correct implant to be used

- VCK
- VCKT
- VCKE
- VCTTi
- VCTTx
- VCTTc
- VCTTcm

- K Implant Visual Card**
- KT Implant Visual Card**
- KE Implant Visual Card**
- TTi Implant Visual Card**
- TTx Implant Visual Card**
- TTc Implant Visual Card**
- TTcm Implant Visual Card**



ST16

**Autoclavable tray**















## Product Code

## Product

### INSTRUMENTS

in surgical steel, sterilisable
















	62SI	Short Screwdriver for internal hex connection Implants
	62MI	Medium Screwdriver for internal hex connection Implants
	62LI	Long Screwdriver for internal hex connection Implants
	99SI	Short Handpiece Screwdriver for internal hex connection Implants
	99MI	Medium Handpiece Screwdriver for internal hex connection Implants
	99LI	Long Handpiece Screwdriver for internal hex connection Implants
	62SX	Short Screwdriver for external hex connection Implants
	62MX	Medium Screwdriver for external hex connection Implants
	62LX	Long Screwdriver for external hex connection Implants
	99SX	Short Handpiece Screwdriver for external hex connection Implants
	99MX	Medium Handpiece Screwdriver for external hex connection Implants
	99LX	Long Handpiece Screwdriver for external hex connection Implants

## Product Code

## Product

### INSTRUMENTS

in surgical steel, sterilisable

	62Sc	Short Screwdriver for conical connection Implants
	62Mc	Medium Screwdriver for conical connection Implants
	62Lc	Long Screwdriver for conical connection Implants
	99Sc	Short Handpiece Screwdriver for conical connection Implants
	99Mc	Medium Handpiece Screwdriver for conical connection Implants
	99Lc	Long Handpiece Screwdriver for conical connection Implants
	62Scm	Short Screwdriver for morse taper connection Implants
	62Mcm	Medium Screwdriver for morse taper connection Implants
	62Lcm	Long Screwdriver for morse taper connection Implants
	99Scm	Short Handpiece Screwdriver for morse taper connection Implants
	99Mcm	Medium Handpiece Screwdriver for morse taper connection Implants
	99Lcm	Long Handpiece Screwdriver for morse taper connection Implants
	99000	Drill Extension Shaft
	60000	Ratchet
	KMP	Adapter for ratchet screwdriver



Product Code

Product

**INSTRUMENTS**

in surgical steel, sterilisable



SDS

Short Manual Screwdriver

SDL

Long Manual Screwdriver

SDTS

Short Torx Manual Screwdriver

SDTL

Long Torx Manual Screwdriver

SDXL

Extra-Long Manual screwdriver



30000S

Short screwdriver with handpiece connection

30000

Screwdriver with handpiece connection

30000TS

Short Torx screwdriver with handpiece connection

30000T

Torx screwdriver with handpiece connection

65000K

**Dynamometric Wrench**

complete with handpiece connection code 66000  
Torque 10-35 N/cm




66000

Handpiece connection

...Protheses  
INTERNAL CONNECTION:  
KAPPA LINE • TORQUE TYPE LINE

# Healing Caps

# Implant Closure Caps



## Product Code

## Product

### TITANIUM HEALING CAPS

#### STRAIGHT

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2
- 5.9 mm. x GH 2
- 6.5 mm. x GH 2

- 3.3 mm. x GH 5
- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5
- 5.9 mm. x GH 5
- 6.5 mm. x GH 5

#### FLARED

- 3.3 mm. x GH 3
- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3
- 5.9 mm. x GH 3
- 6.5 mm. x GH 3

- 3.3 mm. x GH 5
- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5
- 5.9 mm. x GH 5
- 6.5 mm. x GH 5



## Product Code

## Product

### TITANIUM IMPLANT CLOSURE CAPS

- VT33
- VT38
- VT45
- VT52
- VT59

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.



The Healing Cap allows for adequate reconditioning of the peri-implant soft tissue; it should be selected according to the type of prosthetic reconstruction planned.



Product Code

Product

### STANDARD ABUTMENT IMPRESSION COPINGS FOR REPOSITIONING TECHNIQUE

in titanium, cylindrical, VMI prosthetic screw included



IM33  
IM38  
IM45  
IM52  
IM59  
IM65

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.
- 6.5 mm.

### LONG ABUTMENT IMPRESSION COPINGS FOR REPOSITIONING TECHNIQUE

in titanium, flared, VMIL prosthetic screw included




IM33L  
IM38L  
IM45L  
IM52L

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.

### ABUTMENT IMPRESSION COPINGS FOR REPOSITIONING TECHNIQUE

in titanium, VIRB screw included, manual tightening



IRB33  
IRB38  
IRB45  
IRB52

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.



VIRB

Ball-top Screw for short Abutment Impression Copings



VP2

Short Transverse Screw


It is possible to use the VP2 prosthetic screw to replace the VIRB ball-top prosthetic screw.

Product Code

Product

### LONG ABUTMENT IMPRESSION COPINGS FOR PICK-UP TECHNIQUE

in titanium, flared, VMIL prosthetic screw included



IMG33L  
IMG38L  
IMG45L  
IMG52L

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.

## Analogues



A33  
A38  
A45  
A52  
A59  
A65


### Titanium ANALOGUES

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.
- 6.5 mm.


# Prosthetic Accessories

## Synthetic overview

Temporary Prostheses




Cement-retained prosthesis



Single Screw-retained Prosthesis




Multi-unit Screw-retained Prosthesis on EA abutments




directly onto implants

Welding


directly onto implants



CAB on EA abutments



Welding accessories on EA abutments




Overdenture



- I for Internal connection
- X for External connection
- C for Conical connection
- CM for Morse Taper connection







Long connection



Short Connection

SCREW-RETAINED PROSTHESIS  
 BY USING DIFFERENT WINSIX PROSTHETIC COMPONENTS  
 IT IS IMPORTANT TO REMEMBER THAT:

For the screw-retention of prosthetic accessories directly onto the implants  
 only long connection  
 screws may be used  
 (as their length runs through the entire implant hex).



For the use of prosthetic accessories with transmucosal abutments  
 from the FLAT SHIFT line,  
 only short screws may be used.




## Product Code

## Product

### STRAIGHT TEMPORARY ABUTMENTS

in titanium, VP2 prosthetic screw included



- MP3301
- MP3801
- MP4501
- MP5201
- MP5901

- 3.3 mm. x GH 1
- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1
- 5.9 mm. x GH 1



- MP3302
- MP3802
- MP4502
- MP5202
- MP5902

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2
- 5.9 mm. x GH 2

## Product Code

## Product

### WINPeek TEMPORARY ABUTMENTS in peek, with placement indicator on the abutment



WPO332  
WPO382  
WPO452

#### OVAL

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2



WPT332  
WPT382  
WPT452

#### TRIANGULAR

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2



WPF3300  
WPF3800  
WPF4500


#### MILLING

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.

## Product Code

## Product

### STRAIGHT FINAL ABUTMENTS in titanium, VP2 prosthetic screw included



MD3300  
MD3800  
MD4500  
MD5200  
MD5900  
MD6500

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2
- 5.9 mm. x GH 2
- 6.5 mm. x GH 2



MIN3315  
MIN3815  
MIN4515  
MIN5215

### 15° ANGLED FINAL ABUTMENTS in titanium, VP2 prosthetic screw included

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.



MIN3325  
MIN3825  
MIN4525  
MIN5225

### 25° ANGLED FINAL ABUTMENTS in titanium, VP2 prosthetic screw included

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.



MF3300  
MF3800  
MF4500  
MF5200  
MF5900

### MILLING ABUTMENTS in titanium, VP2 prosthetic screw included

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.

## Product Code

## Product

### STRAIGHT MULTIPURPOSE ABUTMENTS

in titanium, VP2 prosthetic screw included



MD3300MP  
MD3800MP  
MD4500MP

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.

### 15° ANGLED MULTIPURPOSE ABUTMENTS

in titanium, VP3 prosthetic screw included




MIN3315MP  
MIN3815MP  
MIN4515MP

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.

### 25° ANGLED MULTIPURPOSE ABUTMENTS

in titanium, VP3 prosthetic screw included



MIN3325MP  
MIN3825MP  
MIN4525MP

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.

## Product Code

## Product

### STRAIGHT ABUTMENTS FOR FINAL TECHNIQUE

in titanium, VP2 prosthetic screw included



MBC3300  
MBC3800  
MBC4500  
MBC5200

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.

### 15° ANGLED ABUTMENTS FOR FINAL TECHNIQUE

in titanium, VP3 prosthetic screw included



MBC3315  
MBC3815  
MBC4515  
MBC5215

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.

### 25° ANGLED ABUTMENTS FOR FINAL TECHNIQUE

in titanium, VP3 prosthetic screw included



MBC3325  
MBC3825  
MBC4525  
MBC5225


- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.

## Product Code

## Product

### STRAIGHT AESTHETIC ABUTMENTS

in titanium, VP2 prosthetic screw included



WE33001  
WE38001  
WE45001

- 3.3 mm. x GH 1 mm.
- 3.8 mm. x GH 1 mm.
- 3.3 mm. x GH 1 mm.

WE33002  
WE38002  
WE45002

- 3.3 mm. x GH 2 mm.
- 3.8 mm. x GH 2 mm.
- 3.3 mm. x GH 2 mm.

WE33003  
WE38003  
WE45003

- 3.3 mm. x GH 3 mm.
- 3.8 mm. x GH 3 mm.
- 3.3 mm. x GH 3 mm.

WE33004  
WE38004  
WE45004


- 3.3 mm. x GH 4 mm.
- 3.8 mm. x GH 4 mm.
- 3.3 mm. x GH 4 mm.

WE33005  
WE38005  
WE45005

- 3.3 mm. x GH 5 mm.
- 3.8 mm. x GH 5 mm.
- 3.3 mm. x GH 5 mm.

### AESTHETIC 15° ANGLED ABUTMENTS

in titanium, VP3 prosthetic screw included



WE33151  
WE38151  
WE45151

- 3.3 mm. x GH 1 mm.
- 3.8 mm. x GH 1 mm.
- 4.5 mm. x GH 1 mm.

WE33152  
WE38152  
WE45152


- 3.3 mm. x GH 2 mm.
- 3.8 mm. x GH 2 mm.
- 4.5 mm. x GH 2 mm.

## Product Code

## Product

### AESTHETIC 15° ANGLED ABUTMENTS

in titanium, VP3 prosthetic screw included



WE33153  
WE38153  
WE45153


- 3.3 mm. x GH 3 mm.
- 3.8 mm. x GH 3 mm.
- 4.5 mm. x GH 3 mm.

WE33154  
WE38154  
WE45154  
WE33155  
WE38155  
WE45155

- 3.3 mm. x GH 4 mm.
- 3.8 mm. x GH 4 mm.
- 4.5 mm. x GH 4 mm.
- 3.3 mm. x GH 5 mm.
- 3.8 mm. x GH 5 mm.
- 4.5 mm. x GH 5 mm.

### AESTHETIC 25° ANGLED ABUTMENTS

in titanium, VP3 prosthetic screw included



WE33251  
WE38251  
WE45251

- 3.3 mm. x GH 1 mm.
- 3.8 mm. x GH 1 mm.
- 4.5 mm. x GH 1 mm.

WE33252  
WE38252  
WE45252

- 3.3 mm. x GH 2 mm.
- 3.8 mm. x GH 2 mm.
- 4.5 mm. x GH 2 mm.

WE33253  
WE38253  
WE45253

- 3.3 mm. x GH 3 mm.
- 3.8 mm. x GH 3 mm.
- 4.5 mm. x GH 3 mm.

WE33254  
WE38254  
WE45254

- 3.3 mm. x GH 4 mm.
- 3.8 mm. x GH 4 mm.
- 4.5 mm. x GH 4 mm.

WE33255  
WE38255  
WE45255


- 3.3 mm. x GH 5 mm.
- 3.8 mm. x GH 5 mm.
- 4.5 mm. x GH 5 mm.



# Prosthetic Abutments

Product Code


Product



MC33  
MC38  
MC45  
MC52  
MC59  
MC65

**CASTABLE ABUTMENTS**  
in PMMA, VP2 prosthetic screw included


- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2
- 5.9 mm. x GH 2
- 6.5 mm. x GH 2



MCF3300  
MCF3800  
MCF4500  
MCF5200  
MCF5900

**CASTABLE DRILLABLE ABUTMENTS**  
in PMMA, VP2 prosthetic screw included

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.



GP33  
GP38  
GP45  
GP52  
GP59

**OVERCASTING ABUTMENTS**  
in PMMA and platinum-gold-palladium-iridium alloys,  
VP2 prosthetic screw included


- 3.3 mm. x GH 1.5
- 3.8 mm. x GH 1.5
- 4.5 mm. x GH 1.5
- 5.2 mm. x GH 1.5
- 5.9 mm. x GH 1.5

...Protheses  
WINZIR and PRO ABUTMENTS  
FOR ZIRCONIA PROSTHESIS

## Product Code

## Product

teardrop shape of the Winzir Abutment seen from above



WZ33  
WZ38  
WZ45  
WZ52  
WZ59

**WINZIR ABUTMENTS**  
in titanium for zirconium crowns,  
VP2 prosthetic screw included


- 3.3 mm. x GH 1
- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1
- 5.9 mm. x GH 1



PRO33  
PRO38  
PRO45

**PRO ABUTMENTS**  
in titanium, for zirconium prostheses,  
VP2 prosthetic screw included

- 3.3 mm. External hex connector
- 3.8 mm. External hex connector
- 4.5 mm. External hex connector



IMPR33  
IMPR38  
IMPR45

**PRO ABUTMENT IMPRESSION COPINGS**  
in titanium, VMIL screw included

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.



VP6

**Long Prosthetic Screw for PRO**



use with VP6

# ...Prostheses

## OVERDENTURE COMPONENTS

# LOCATOR Attachments

for the correct and durable placement of Overdentures



## for Internal Hex Implants

02701	Ø 3.3 - H 1 mm		02712	Ø 4.5 - H 1 mm	
02702	Ø 3.3 - H 2 mm		02713	Ø 4.5 - H 2 mm	
02703	Ø 3.3 - H 3 mm		02714	Ø 4.5 - H 3 mm	
02704	Ø 3.3 - H 4 mm		02715	Ø 4.5 - H 4 mm	
02705	Ø 3.3 - H 5 mm		02716	Ø 4.5 - H 5 mm	
02707	Ø 3.8 - H 1 mm		02718	Ø 5.2 - H 1 mm	
02708	Ø 3.8 - H 2 mm		02719	Ø 5.2 - H 2 mm	
02709	Ø 3.8 - H 3 mm		02720	Ø 5.2 - H 3 mm	
02710	Ø 3.8 - H 4 mm		02721	Ø 5.2 - H 4 mm	
02711	Ø 3.8 - H 5 mm		02722	Ø 5.2 - H 5 mm	


## Instrument Product List

Product Code	Product	
08393	LOCATOR Core Tool	
08913	LOCATOR driver with CA connection - 23 mm	
08914	LOCATOR driver with CA connection - 29 mm	
08517	Parallel Post (4 Pcs)	
09530	Angle measurement guide	
08505	LOCATOR Abutment Impression Coping (4 Pcs)	
08530	LOCATOR Analogue Ø 4.0 mm (4 Pcs)	
08516	LOCATOR Analogue Ø 5.0 mm (4 Pcs)	

## Versatility of Retention and Angulation


### Retainers for disparallelism up to 20°

For prosthetic application on angled implants up to 10°, allowing a divergence between two implants of up to 20°.



### Retainers for disparallelism up to 40°




For prosthetic application on angled implants up to 20°, allowing a divergence between two implants of up to 40°.



08515	Black processing replacement male	
08510-V	Titanium metal cap (4 Pcs)	
08514	Block-out spacer - silicon rubber (white)	
08519-2	Assorted retainers with titanium cap, disparallelism up to 20° (2 Pcs each)	
08550-2	Assorted retainers with steel cap, disparallelism up to 20° (2 Pcs each)	
08524	Transparent retainers, disparallelism up to 20° - 22.25 N (4 Pcs)	
08527	Pink Retainers, disparallelism up to 20° - 13.36 N (4 Pcs)	
08529	Blue retainers, disparallelism up to 20° - 6.68 N (4 Pcs)	
08540-2	Assorted retainers with titanium cap, disparallelism up to 40° (2 Pcs each)	
08547	Green Retainers, disparallelism up to 40° - 17.80 N (4 Pcs)	
08548	Red Retainers, disparallelism up to 40° - 4.45 N (4 Pcs)	
08915	Orange Retainers, disparallelism up to 40° - 8.90 N (4 Pcs)	






PRODUCT FROM: ZEST ANCHORS

# Ball Attachments

Product Code	Product
 MOR331 MOR381 MOR451	<b>BALL ATTACHMENTS</b> consisting of: O-ring Abutment and titanium Box, O-ring in NBR, ball Ø 2.1 mm. <ul style="list-style-type: none"> <li>● 3.3 mm. x GH 1</li> <li>● 3.8 mm. x GH 1</li> <li>● 4.5 mm. x GH 1</li> </ul>
 MOR333 MOR383 MOR453	<ul style="list-style-type: none"> <li>● 3.3 mm. x GH 3</li> <li>● 3.8 mm. x GH 3</li> <li>● 4.5 mm. x GH 3</li> </ul>
 KR 99KR	<b>WRENCHES</b> <b>Ball attachment wrench</b> <b>Wrench for ball attachments with handpiece connection</b>

N.B. To secure MOR Ball attachments a KR prosthetic wrench is required, to be purchased separately.

## Accessories for Ball Attachments

	BOR	Box for ball attachment complete with OR
	OR	O-ring
	ORK	Special O-ring
<b>MICRO BOR</b>		
	Micro BOR	Box for ball attachment complete with micro OR
	Micro OR	Micro O-ring



# ...Prostheses

FLAT SHIFT LINE EA EXTREME  
 ABUTMENT LINE MULTI-UNIT  
 SCREW-RETAINED PROSTHESIS



In the case of MULTI-UNIT SCREW-RETAINED PROSTHESES, with implant disparallelism of up to 45°, the WINSIX System provides the Flat Shift Prosthetic Line, which allows the prosthesis to be placed directly onto the implant, as well as onto the Flat Shift transmucosal abutment in the event of very dense soft tissue. The Flat Shift Line also allows the construction of direct and passive screw-retained prostheses even in cases of strong disparallelism.

## IMPRESSION COPINGS AND PROSTHETIC ABUTMENTS




Product Code	Product
	<b>FLAT SHIFT ABUTMENTS</b> in titanium, VP2 prosthetic screw included
MFS333 MFS383 MFS453 MFS523	<ul style="list-style-type: none"> <li>● 3.3 mm. x GH 3</li> <li>● 3.8 mm. x GH 3</li> <li>● 4.5 mm. x GH 3</li> <li>● 5.2 mm. x GH 3</li> </ul>
MFS334 MFS384 MFS454 MFS524	<ul style="list-style-type: none"> <li>● 3.3 mm. x GH 4</li> <li>● 3.8 mm. x GH 4</li> <li>● 4.5 mm. x GH 4</li> <li>● 5.2 mm. x GH 4</li> </ul>
MFS335 MFS385 MFS455 MFS525	<ul style="list-style-type: none"> <li>● 3.3 mm. x GH 5</li> <li>● 3.8 mm. x GH 5</li> <li>● 4.5 mm. x GH 5</li> <li>● 5.2 mm. x GH 5</li> </ul>
	<b>FLAT SHIFT CLOSURE CAPS</b> in titanium
VTFS33 VTFS38 VTFS45 VTFS52	<ul style="list-style-type: none"> <li>● 3.3 mm.</li> <li>● 3.8 mm.</li> <li>● 4.5 mm.</li> <li>● 5.2 mm.</li> </ul>
VFS	<b>Flat Shift Prosthetic screws</b>

## Abutment Impression Copings

Product Code


Product



IFS33  
IFS38  
IFS45  
IFS52

**FLAT SHIFT ABUTMENT IMPRESSION COPINGS FOR REPOSITIONING TECHNIQUE**  
in titanium, VIFS screw included

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.



IFS33L  
IFS38L  
IFS45L

**FLAT SHIFT ABUTMENT IMPRESSION COPINGS FOR PICK-UP TECHNIQUE**  
in titanium, flared, VMIL screw included

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.

Product Code

Product



VIFS

**ABUTMENT SCREWS**

**Screw for Flat Shift Abutment Impression Coping**

The VIFS Screw cannot be used directly with the implant, but can only be applied to the Flat Shift Abutment



VMI

**Screw for Short Abutment Impression Coping**



VMIL

**Screw for Long Abutment Impression Coping**

To take the impression of the implant directly with the Abutment Impression Coping Flat Shift for the repositioning technique requires the use of VMI or VMIL Screws

## Product Code

## Product

### RETENTIVE FLARED DIRECT ABUTMENTS FLAT SHIFT

in titanium, VP2 prosthetic screw included




MDS331WR  
MDS381WR  
MDS451WR  
MDS521WR

- 3.3 mm. x GH 1
- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1



The MDS Abutments can also be used for welding either directly in the mouth of the patient, or possibly in the laboratory in the case of temporary immediate loading prostheses



MDS332WR  
MDS382WR  
MDS452WR  
MDS522WR

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2



MDS333WR  
MDS383WR  
MDS453WR  
MDS523WR

- 3.3 mm. x GH 3
- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3

## Product Code

## Product

### CYLINDRICAL DIRECT ABUTMENTS FLAT SHIFT

in titanium, VP2 prosthetic screw included




MDC331  
MDC381  
MDC451  
MDC521

- 3.3 mm. x GH 1
- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1



MDC332  
MDC382  
MDC452  
MDC522

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2



MDC333  
MDC383  
MDC453  
MDC523

- 3.3 mm. x GH 3
- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3

# Flat Shift Line

## Product Code

## Product

### FLAT SHIFT CASTABLE CAPS

In PMMA, for cylindrical direct abutment




CMD33F  
CMD38F  
CMD45F  
CMD52F

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.

### CASTABLE DIRECT ABUTMENTS FLAT SHIFT

in PMMA, VP2 prosthetic screw included



MDF331  
MDF381  
MDF451  
MDF521


- 3.3 mm. x GH 1
- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1



# EA LINE - Extreme Abutment

The EA - Extreme Abutment - Line of prosthetic accessories has been specially designed for those cases involving a multi-unit screw-retained prosthetic restoration. They allow the placement of WINSIX implants in a tilted position for a restoration on permanent and stable implants.

EAAbutments, both straight and with angulation of 20° and 30° - find their ideal application in the innovative Just on 4/6 Techniques, enabling the prosthetic application of implants even with a strong inclination of up to 45°.



## IMPRESSION COPINGS AND PROSTHETIC ABUTMENTS

### EA LINE - Extreme Abutment

WINSIX Implants

EAT  
Abutment Impression Copings for EA

EA straight      EA angled

EAS  
Smooth  
Abutment

EAF  
Castable  
Abutment

EAR  
Retentive  
Abutment

EAA  
Analogue



Preassembled Mouter included in the EA straight abutment package for correct placement and initial tightening laps.



Tightening of the Abutment must then be carried out at 30 N/cm with the Handpiece Screwdriver code 99EA.

Manual screwdriver code 62EA.



internal

## Prosthetic Abutments

Product Code

Product

**EA straight ABUTMENTS**  
in titanium,  
Preassembled Mouter included



1 mm.

EA33001  
EA38001  
EA45001  
EA52001

- 3.3 mm. x GH 1
- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1

**EA straight ABUTMENTS**  
in titanium,  
Preassembled Mouter included



2 mm.

EA33002  
EA38002  
EA45002  
EA52002

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2

# EA LINE - Extreme Abutment



Preassembled Mounter included in the EA straight abutment package for correct placement and initial tightening laps.



Tightening of the Abutment must then be carried out at 30 N/cm with the Handpiece Screwdriver code 99EA.

Manual screwdriver code 62EA.




internal

## Product Code

## Product


**EA straight ABUTMENTS**  
in titanium,  
Preassembled Mounter included



EA33003  
EA38003  
EA45003  
EA52003

- 3.3 mm. x GH 3
- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3

**EA straight ABUTMENTS**  
in titanium,  
Preassembled Mounter included




EA33004  
EA38004  
EA45004  
EA52004

- 3.3 mm. x GH 4
- 3.8 mm. x GH 4
- 4.5 mm. x GH 4
- 5.2 mm. x GH 4

## Product Code

## Product

**EA straight ABUTMENTS**  
in titanium,  
Preassembled Mounter included



EA33005  
EA38005  
EA45005  
EA52005

- 3.3 mm. x GH 5
- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5

## Product Code

## Product


### EA 20° ANGLED ABUTMENTS

in titanium,  
VP3 prosthetic screw included,  
Preassembled Mounter included




EA33202  
EA38202  
EA45202  
EA52202

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2



EA33203  
EA38203  
EA45203  
EA52203

- 3.3 mm. x GH 3
- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3



EA33204  
EA38204  
EA45204  
EA52204

- 3.3 mm. x GH 4
- 3.8 mm. x GH 4
- 4.5 mm. x GH 4
- 5.2 mm. x GH 4



EA Angled Mounter

**Preassembled mounter**  
included in the angled  
Abutments package for placement  
and verification of the angulation.

## Product Code

## Product

### EA 30° ANGLED ABUTMENTS

in titanium,  
VP3 prosthetic screw included,  
Preassembled Mounter included



EA33304  
EA38304  
EA45304  
EA52304

- 3.3 mm. x GH 4
- 3.8 mm. x GH 4
- 4.5 mm. x GH 4
- 5.2 mm. x GH 4



EA33305  
EA38305  
EA45305  
EA52305

- 3.3 mm. x GH 5
- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5



EA Angled Mounter

**Preassembled mounter**  
included in the angled  
Abutments package for placement  
and verification of the angulation.

Product Code

Product

**EA ACCESSORIES AND PROSTHETIC INSTRUMENTS**



EAA

**Analogue for EA**  
in titanium



EAT

**Abutment Impression Coping for EA**  
in titanium, VEAT prosthetic screw included



EAR

**Retentive Abutment for EA**  
in titanium, VP4 prosthetic screw included



EAS

**Smooth Abutment for EA**  
in titanium, VP4 prosthetic screw included



EAF

**EAF Castable Abutment for EA**  
in PMMA, VP4 prosthetic screw included

Product Code

Product

**EA ACCESSORIES AND PROSTHETIC INSTRUMENTS**



VGEA

**Healing Cap for EA**  
in titanium



VP3

**Screw for EA Connection - Implant**



VP4

**Prosthetic screw for EA/EAx Abutment**



VEAT

**Screw for Abutment Impression Coping for EA**



99EA

**Straight EA screwdriver for handpiece**













62EA

**Straight EA manual screwdriver**



# WINSIX Screws

internal hex

Product Code	Product
	<b>WINSIX PROSTHESIC SCREWS</b> in titanium
 VP1	<b>Long Transverse screw</b>
 VP2	<b>Short Transverse Screw</b>
 VP3	<b>VP3 Screw for EA connection</b>
 VP4	<b>Prosthetic screw for EA Abutment</b>
 VP6	<b>Long Prosthetic Screw for PRO</b>
 VMI	<b>Screw for Short Abutment Impression Coping</b>
 VMIL	<b>Screw for Long Abutment Impression Coping</b>
 VIRB	<b>Screw for short Abutment Impression Coping</b>
 VFS	<b>Prosthetic screw for Flat Shift Abutment</b>
 VIFS	<b>Screw for Flat Shift Abutment Impression Coping</b>
 VEAT	<b>Screw for Abutment Impression Coping for EA</b>

TORQUE VALUES TO BE APPLIED DURING POSITIONING OF THE SCREWS USING THE DYNAMOMETRIC WRENCH	
VT e VG	⊙ 10 N/cm
VP4	⊙ 20 N/cm
VFS	⊙ 25 N/cm
VP2	⊙ 30 N/cm
VP3	⊙ 30 N/cm
EA straight	⊙ 30 N/cm
MFS ABUTMENT	⊙ 30 N/cm

...Protheses  
EXTERNAL CONNECTION:  
TTx TORQUE TYPE

## Product Code

## Product

### HEALING CAPS

in titanium, to be sterilised before use



VG3303x  
VG3803x  
VG4503x  
VG5203x  
VG5903x

#### STRAIGHT

- 3.3 mm. x GH 3
- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3
- 5.9 mm. x GH 3

The straight healing caps are 5 mm in height they can also be used as a guide for the Bone mill BM series. (See page 11)



VG3305x  
VG3805x  
VG4505x  
VG5205x  
VG5905x

- 3.3 mm. x GH 5
- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5
- 5.9 mm. x GH 5



VG333x  
VG383x  
VG453x  
VG523x  
VG593x

#### FLARED

- 3.3 mm. x GH 3
- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3
- 5.9 mm. x GH 3

VG335x  
VG385x  
VG455x  
VG525x  
VG595x

- 3.3 mm. x GH 5
- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5
- 5.9 mm. x GH 5

### IMPLANT CLOSURE CAPS

in titanium, to be sterilised before use



VT33x  
VT38x  
VT45x  
VT52x  
VT59x

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.

## Product Code

## Product

### ABUTMENT IMPRESSION COPINGS FOR PICK-UP TECHNIQUE


with VMILx prosthetic screw included



IMG33Lx  
IMG38Lx  
IMG45Lx  
IMG52Lx  
IMG59Lx

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.

## Analogues



A33x  
A38x  
A45x  
A52x  
A59x

### ANALOGUES

in titanium

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.

## Product Code

## Product

### STRAIGHT TEMPORARY ABUTMENTS

in titanium, VP2x prosthetic screw included



MP3301x  
MP3801x  
MP4501x  
MP5201x  
MP5901x

- 3.3 mm. x GH 1
- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1
- 5.9 mm. x GH 1




MP3302x  
MP3802x  
MP4502x  
MP5202x  
MP5902x

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2
- 5.9 mm. x GH 2

### ROTATIONAL TEMPORARY ABUTMENTS

in titanium, VP2x prosthetic screw included



MP3301RX  
MP3801RX

- 3.3 mm. x GH 1
- 3.8 mm. x GH 1

### CYLINDRICAL DIRECT ABUTMENTS

in titanium, VP2x prosthetic screw included



MDC331X  
MDC381X


- 3.3 mm. x GH 1
- 3.8 mm. x GH 1

## Product Code

## Product

### FINAL ABUTMENTS

in titanium, VP2x prosthetic screw included




MD3300x  
MD3800x  
MD4500x  
MD5200x  
MD5900x

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2
- 5.9 mm. x GH 2

### 15° ANGLED FINAL ABUTMENTS

in titanium, VP2x prosthetic screw included



MIN3315x  
MIN3815x  
MIN4515x  
MIN5215x  
MIN5915x

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.

### 25° ANGLED FINAL ABUTMENTS

in titanium, VP2x prosthetic screw included



MIN3325x  
MIN3825x  
MIN4525x  
MIN5225x  
MIN5925x

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.

### MILLING ABUTMENTS

in titanium, VP2x prosthetic screw included



MF3300x  
MF3800x  
MF4500x  
MF5200x  
MF5900x

- 3.3 mm.
- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.


# Prosthetic Abutments

Product Code

Product

## CASTABLE ABUTMENTS

in PMMA, VP2x prosthetic screw included




MC33x  
MC38x  
MC45x  
MC52x  
MC59x

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2
- 5.9 mm. x GH 2

## OVERCASTING ABUTMENTS

IN PMMA AND PLATINOR  
VP2x prosthetic screw included



GP33x  
GP38x  
GP45x  
GP52x  
GP59x

- 3.3 mm. x GH 1.5
- 3.8 mm. x GH 1.5
- 4.5 mm. x GH 1.5
- 5.2 mm. x GH 1.5
- 5.9 mm. x GH 1.5

## CASTABLE ABUTMENTS

in PMMA, VP2x prosthetic screw included



C33x  
C38x  
C45x


- 3.3 mm.
- 3.8 mm.
- 4.5 mm.

...Protheses  
MULTI-UNIT SCREW-RETAINED  
EAX LINE - EXTREME ABUTMENT



# EAX Line - Extreme Abutment

EA Abutments, both straight and with angulation of 17° and 30° - find their ideal application in the innovative Just on 4/6 Techniques, enabling the prosthetic application of implants even with a strong inclination of up to 45°.



Preassembled Mouter included in the EA straight abutment package for correct placement and initial tightening laps.







Tightening of the Abutment must then be carried out at 30 N/cm with the Handpiece Screwdriver code 99EA

Manual screwdriver code 62EA




## EAX Prosthetic Abutments

Product Code	Product
 EA33001x EA38001x EA45001x EA52001x	<b>EAX straight ABUTMENTS</b> in titanium, Preassembled Mouter included  ● 3.3 mm. x GH 1 ● 3.8 mm. x GH 1 ● 4.5 mm. x GH 1 ● 5.2 mm. x GH 1
 EA33002x EA38002x EA45002x EA52002x	<b>EAX straight ABUTMENTS</b> in titanium, Preassembled Mouter included  ● 3.3 mm. x GH 2 ● 3.8 mm. x GH 2 ● 4.5 mm. x GH 2 ● 5.2 mm. x GH 2
 EA33003x EA38003x EA45003x EA52003x	<b>EAX straight ABUTMENTS</b> in titanium, Preassembled Mouter included  ● 3.3 mm. x GH 3 ● 3.8 mm. x GH 3 ● 4.5 mm. x GH 3 ● 5.2 mm. x GH 3



EA33004x  
EA38004x  
EA45004x  
EA52004x

- 3.3 mm. x GH 4
- 3.8 mm. x GH 4
- 4.5 mm. x GH 4
- 5.2 mm. x GH 4



EA33005x  
EA38005x  
EA45005x  
EA52005x

**EAX straight ABUTMENTS**  
 in titanium,  
 Preassembled Mouter included

- 3.3 mm. x GH 5
- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5

## Product Code

## Product

### EAx 17° ANGLED ABUTMENTS

in titanium, VP3x prosthetic screw included  
Preassembled Mounter included



EA33172x  
EA38172x  
EA45172x  
EA52172x

- 3.3 mm. x GH 2
- 3.8 mm. x GH 2
- 4.5 mm. x GH 2
- 5.2 mm. x GH 2

### EAx 17° ANGLED ABUTMENTS

in titanium, VP3x prosthetic screw included  
Preassembled Mounter included



EA33173x  
EA38173x  
EA45173x  
EA52173x

- 3.3 mm. x GH 3
- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3

### EAx 17° ANGLED ABUTMENTS

in titanium, VP3x prosthetic screw included  
Preassembled Mounter included



EA33174x  
EA38174x  
EA45174x  
EA52174x

- 3.3 mm. x GH 4
- 3.8 mm. x GH 4
- 4.5 mm. x GH 4
- 5.2 mm. x GH 4



EA Angled Mounter

**Preassembled mounter**  
included in the angled Abutments  
package for placement and  
verification of the angulation.

## Product Code

## Product

### EAX 30° ANGLED ABUTMENTS

in titanium, VP3x prosthetic screw included  
Preassembled Mounter included



EA33304x  
EA38304x  
EA45304x  
EA52304x

- 3.3 mm. x GH 4
- 3.8 mm. x GH 4
- 4.5 mm. x GH 4
- 5.2 mm. x GH 4

### EAX 30° ANGLED ABUTMENTS °

in titanium, VP3x prosthetic screw included  
Preassembled Mounter included



EA33305x  
EA38305x  
EA45305x  
EA52305x







- 3.3 mm. x GH 5
- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5








EA Angled Mounter

**Preassembled mounter**  
included in the angled Abutments  
package for placement and  
verification of the angulation.



Product Code	Product
	<b>EAx PROSTHETIC ACCESSORIES</b> in titanium
 EAAx	<b>Analogue for EAx</b> in titanium
 EATx	<b>Abutment Impression Coping for EAx</b> in titanium, VEATx impression screws included
 12 mm. EARx	<b>Retentive Abutment for EAx</b> in titanium, VP4 prosthetic screw included
 10 mm. EASx	<b>Smooth Abutment for EAx</b> in titanium, VP4 prosthetic screw included
 8 mm. EAFx	<b>Castable Abutment for EAx</b> in PMMA, VP4 prosthetic screw included
 VGEAx	<b>Healing Cap for EAx</b> in titanium

Product Code	Product
	<b>WINSIX PROSTHETIC SCREWS</b> in titanium
 VP2x	<b>Transverse Screw</b>
 VP3x	<b>Screw for EAx Connection - Implant</b>
 VP4	<b>Abutment Screw for EA/EAx</b>
 VMILx	<b>Screw for Long Abutment Impression Coping</b>
 VEATx	<b>Screw for Abutment Impression Coping for EAx</b>

TORQUE VALUES TO BE APPLIED DURING POSITIONING OF THE SCREWS USING THE DYNAMOMETRIC WRENCH	
VT e VG	C 10 N/cm
VP4	C 20 N/cm
VP2X	C 30 N/cm
VP3X	C 30 N/cm
EA straight	C 30 N/cm

# ...Prostheses

OVERDENTURE COMPONENTS

# LOCATOR Attachments

for the correct and durable placement of Overdentures



## For External Hex Implants

02865	Ø 3.3 - H 1 mm		02877	Ø 4.5 - H 1 mm	
02866	Ø 3.3 - H 2 mm		02878	Ø 4.5 - H 2 mm	
02867	Ø 3.3 - H 3 mm		02879	Ø 4.5 - H 3 mm	
02868	Ø 3.3 - H 4 mm		02880	Ø 4.5 - H 4 mm	
02869	Ø 3.3 - H 5 mm		02881	Ø 4.5 - H 5 mm	
02871	Ø 3.8 - H 1 mm		02883	Ø 5.2 - H 1 mm	
02872	Ø 3.8 - H 2 mm		02884	Ø 5.2 - H 2 mm	
02873	Ø 3.8 - H 3 mm		02885	Ø 5.2 - H 3 mm	
02874	Ø 3.8 - H 4 mm		02886	Ø 5.2 - H 4 mm	
02875	Ø 3.8 - H 5 mm		02887	Ø 5.2 - H 5 mm	


## Instrument Product List

Product Code	Product	
08393	LOCATOR Core Tool	
08913	LOCATOR driver with CA connection - 23 mm	
08914	LOCATOR driver with CA connection - 29 mm	
08517	Parallel Post (4 Pcs)	
09530	Angle measurement guide	
08505	LOCATOR Abutment Impression Coping (4 Pcs)	
08530	LOCATOR Analogue Ø 4.0 mm (4 Pcs)	
08516	LOCATOR Analogue Ø 5.0 mm (4 Pcs)	

## Versatility of Retention and Angulation


### Retainers for disparallelism up to 20°

For prosthetic application on angled implants up to 10°, allowing a divergence between two implants of up to 20°.



### Retainers for disparallelism up to 40°

For prosthetic application on angled implants up to 20°, allowing a divergence between two implants of up to 40°.



08515	Black processing replacement male	
08510-V	Titanium metal cap (4 Pcs)	
08514	Block-out spacer - silicon rubber (white)	
08519-2	Assorted retainers with titanium cap, disparallelism up to 20° (2 Pcs each)	
08550-2	Assorted retainers with steel cap, disparallelism up to 20° (2 Pcs each)	
08524	Transparent retainers, disparallelism up to 20° - 22.25 N (4 Pcs)	
08527	Pink Retainers, disparallelism up to 20° - 13.36 N (4 Pcs)	
08529	Blue retainers, disparallelism up to 20° - 6.68 N (4 Pcs)	
08540-2	Assorted retainers with titanium cap, disparallelism up to 40° (2 Pcs each)	
08547	Green Retainers, disparallelism up to 40° - 17.80 N (4 Pcs)	
08548	Red Retainers, disparallelism up to 40° - 4.45 N (4 Pcs)	
08915	Orange Retainers, disparallelism up to 40° - 8.90 N (4 Pcs)	

PRODUCT FROM: ZEST ANCHORS



# ...Protheses

CONICAL CONNECTION  
TTC IMPLANTS

## Product Code

## Product

### HEALING CAPS in titanium



VG381C  
VG451C  
VG521C

#### FLARED

- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1



VG383C  
VG453C  
VG523C

- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3



VG385C  
VG455C  
VG525C

- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5

### IMPLANT CLOSURE CAPS in titanium



VT38C  
VT45C  
VT52C  
VT59C

- 3.8 mm.
- 4.5 mm.
- 5.2 mm.
- 5.9 mm.

## Product Code

## Product

### STANDARD ABUTMENTS FOR REPOSITIONING TECHNIQUE in titanium, VMICc screw included



IM38c

- 3.8 mm.



IMG38Lc

### LONG ABUTMENTS FOR PICK-UP TECHNIQUE in titanium, VMILc screw included

- 3.8 mm.

## Analogues



A38C

### ANALOGUES in titanium

- 3.8 mm.

## Product Code

## Product

### STRAIGHT TEMPORARY ABUTMENTS

in titanium, VP2c prosthetic screw included



MP3803C  
MP4503C  
MP5203C

- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3



MP3805C  
MP4505C  
MP5205C

- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5

## Product Code

## Product

### STRAIGHT FINAL ABUTMENTS

in titanium, VP2c prosthetic screw included



MD3801C  
MD4501C  
MD5201C

- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1

MD3803C  
MD4503C  
MD5203C

- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3

MD3805C  
MD4505C  
MD5205C

- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5

### 15° ANGLED FINAL ABUTMENTS

in titanium, VP2c prosthetic screw included



MIN38151C  
MIN45151C  
MIN52151C

- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1

MIN38153C  
MIN45153C  
MIN52153C

- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3

MIN38155C  
MIN45155C  
MIN52155C

- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5



## Product Code

## Product

### 25° ANGLED FINAL ABUTMENTS

in titanium, VP2c prosthetic screw included



MIN38251C  
MIN45251C  
MIN52251C

- 3.8 mm. x GH 1
- 4.5 mm. x GH 1
- 5.2 mm. x GH 1

MIN38253C  
MIN45253C  
MIN52253C

- 3.8 mm. x GH 3
- 4.5 mm. x GH 3
- 5.2 mm. x GH 3

MIN38255C  
MIN45255C  
MIN52255C

- 3.8 mm. x GH 5
- 4.5 mm. x GH 5
- 5.2 mm. x GH 5

### DRILLABLE ABUTMENTS

in titanium, VP2c prosthetic screw included



MF3800C  
MF4500C  
MF5200C

- 3.8 mm.
- 4.5 mm.
- 5.2 mm.

## Product Code

## Product

### MULTIFUNCTIONAL STRAIGHT

in titanium, prosthetic VTMF and VPGc screw included



EA003MFC  
EA004MFC  
EA005MFC

GH 3  
GH 4  
GH 5

### MULTIFUNCTIONAL ANGLED 17°

in titanium, VTMF and VPGc prosthetic screw included



EA173MFC  
EA174MFC  
EA175MFC

GH 3  
GH 4  
GH 5

### MULTIFUNCTIONAL ANGLED 30°

in titanium, VTMF and VPGc prosthetic screw included



EA303MFC  
EA304MFC  
EA305MFC

GH 3  
GH 4  
GH 5

VPGC

### VPGC PROSTHETIC SCREW

in titanium for extreme Multifunctional abutment



VTMF

### CLOSURE CAP

in titanium, for extreme Multifunctional Abutment



# Multifunctional Abutments



EACONE

**MULTIFUNCTIONAL**  
in titanium



EAW

**CONICAL ATTACHMENT**  
titanium, for extreme Multifunctional abutment.  
Can be used with: VGEAx - EATx - EARx - EASx - EAFx - AT



EACONEW

**CONICAL ATTACHMENT FOR EAW**  
in titanium, for extreme Multifunctional abutment



EAMOR

**BALL ATTACHMENT**  
in titanium, for extreme Multifunctional Abutment  
with titanium Box and O-Ring

## Through Screws



VMICc

**Screw for Short Abutment Impression Coping**

VMILc

**Screw for Long Abutment Impression Coping**

VP2c

**Short transverse screw**

VPGC

**Multifunctional prosthetic screw**

VTMF

**Multifunctional closure cap**

TORQUE VALUES TO BE APPLIED DURING POSITIONING OF THE SCREWS USING THE DYNAMOMETRIC WRENCH	
VTc	C 10 N/cm
VTMF	C 10 N/cm
VGc	C 10 N/cm
VP2c	C 20 N/cm
VPGc	C 20 N/cm
VP4	C 20 N/cm

# ...Prostheses

CAB DEVICE

CLIP ABUTMENT BAR





# New Restorative Techniques for Multi-Unit Screw-Retained Prostheses

## CAB Clip Abutment Bar



### Teeth Just on 4

In order to offer prosthetically satisfactory and economically viable restorative solutions to a significant portion of the population, WINSIX has developed surgical and prosthetic protocols for resolving cases for patients who are severely edentulous in the upper or lower jaws:



**Teeth Just on 4th and Teeth Just on 6.**

Clinical studies and research have shown that the inclined placement of two distal implants is an effective and simple technique to compensate for any bone failure, thus expanding the prosthetic base in a stable and functional manner. In the case of the mandible, this technique also protects final portion of the arch, avoiding interference with the mandibular nerve.

The **Teeth Just On 4 or Teeth Just on 6** technique permits complete, permanent and stable arch restoration, in many cases avoiding the necessity of bone regeneration interventions and the resulting discomfort and cost to the patient involved.

The new **CAB CLIP ABUTMENT BAR** device permits solidarization of dental implants by means of a Clip - Abutment - Bar combination applied to immediate and non-immediate loading, screw-retained prostheses, in line with the Teeth Just on 4/6 restorative technique protocol.

The CAB is used to build a passive titanium structure, to be implemented in the shortest time, in the context of immediate load application. It represents a reinforcement of the temporary prosthesis, minimising the risk of fracture that may lead to failure of the Implants.



Teeth just on **4**

**MANDIBLE: TEETH JUST ON 4**  
For the complete, permanent and stable prosthetic restoration of the lower arch on only 4 implants.



The patent issued for the CAB® by the United States authority confirms its originality and innovative technology content.  
 EP Patent no. 11425032.7  
 PCT/EP2011/072448  
 EP Patent n.11425032.7

# CAB Clip Abutment Bar



## Bar preparation procedure

After the placement of the EAx Line implants and their respective abutments, the impression is taken with the appropriate EATx abutment impression copings. Following this, the laboratory model with analogues for EAx placement is developed, the AT abutments are selected according to the shoulder height, allowing the bar to be parallel to the occlusal plane. Once the ATxx abutments are placed, each individual bar is cut using a dedicated CB instrument.

After cutting all the bars, the dedicated CF1 or CM1 Clip is inserted, and the assembled structure is mounted on the ATxx abutments and permanently fixed to them by means of specific cements.

## Kit CAB



## CAB Components

<b>Code BTO</b> Bar without Clip	
<b>Code BT1</b> Retentive bar with Male clip, cylindrical hole and 1 mm thickness	
<b>Code BT2</b> Retentive bar equipped with male clip with cylindrical hole and 2 mm thickness	
<b>Code BT3</b> Retentive bar equipped with male clip with elliptical hole and 2 mm thickness	
<b>Code CM1</b> Male Clip with cylindrical hole and 1 mm thickness	
<b>Code CM2</b> Male Clip with cylindrical hole and 2 mm thickness	
<b>Code CM3</b> Male Clip with elliptical hole and 2 mm thickness	
<b>Code CF1</b> Female Clip	
<b>Code AT.</b> Variable abutment shoulder where xx indicates variable shoulder height from 1.7 to 4.2 mm	
<b>Code CB</b> Cutter Bar	




The CAB can be used with the most widely used implant systems with external hex connection.

## Fully adjustable bar

Thanks to the cylindrical or elliptical geometry of the special clips, the Clip-Abutment connection allows extreme versatility of use even in the event of severe disparallelism. This feature makes it ideal for multi-unit screw-retained prostheses with immediate or deferred loading.



1. Measurement of the CAB bar with the help of the Cutter Bar device



2. Cutting the CAB bar with separating disc



3. The mounted CAB bar is easily positioned, even in the event of severe disparallelism

# CAB Surgical kit configuration

To carry out the first CAB case you must have a complete CAB Kit.  
 After the first case, there will be a minimum of 20 to a maximum of 50 Abutments available in the KIT.  
 In the following cases, it will be sufficient to reorder only the missing equipment.  
 All kits include ATxx Abutments for each shoulder height from 1.7 mm to 4.2 mm.

Just On 4 KIT		
Code	Description	Quantity
CM1	Male Clip - Cylindrical Hole 1 mm Thickness	1
CF1	Female Clip	2
BT1	Bar + Male Clip Cylindrical Hole 1 mm thickness	1
BT2	Bar + Male Clip Cylindrical Hole 2 mm thickness	2
BT3	Bar + Male Clip Elliptical Hole 2 mm thickness	2
AT17	Abutment Shoulder 1.7 mm	4
AT22	Abutment Shoulder 2.2 mm	4
AT27	Abutment Shoulder 2.7 mm	4
AT32	Abutment Shoulder 3.2 mm	4
AT37	Abutment Shoulder 3.7 mm	4
AT42	Abutment Shoulder 4.2 mm	4
EAAx	EAx Analogue line	4
KMP	Manual Adapter for Ratchet Instruments	2
62T	Torx wrench with ratchet attachment	2
CB	Cutter Bar	2
Tray	Instrument Caddy	1

Just On 6 KIT		
Code	Description	Quantity
CM1	Male Clip - Cylindrical Hole 1 mm Thickness	1
CF1	Female Clip	4
BT1	Bar + Male Clip Cylindrical Hole 1 mm thickness	3
BT2	Bar + Male Clip Cylindrical Hole 2 mm thickness	2
BT3	Bar + Male Clip Elliptical Hole 2 mm thickness	2
AT17	Abutment Shoulder 1.7 mm	6
AT22	Abutment Shoulder 2.2 mm	6
AT27	Abutment Shoulder 2.7 mm	6
AT32	Abutment Shoulder 3.2 mm	6
AT37	Abutment Shoulder 3.7 mm	6
AT42	Abutment Shoulder 4.2 mm	6
EAAx	EAx Analogue line	6
KMP	Manual Adapter for Ratchet Instruments	2
62T	Torx wrench with ratchet attachment	2
CB	Cutter Bar	2
Tray	Instrument Caddy	1

JUST ON 4/6 KIT		
Code	Description	Quantity
CM1	Male Clip - Cylindrical Hole 1 mm Thickness	2
CF1	Female Clip	6
BT1	Bar + Male Clip Cylindrical Hole 1 mm thickness	4
BT2	Bar + Male Clip Cylindrical Hole 2 mm thickness	4
BT3	Bar + Male Clip Elliptical Hole 2 mm thickness	4
AT17	Abutment Shoulder 1.7 mm	10
AT22	Abutment Shoulder 2.2 mm	10
AT27	Abutment Shoulder 2.7 mm	10
AT32	Abutment Shoulder 3.2 mm	10
AT37	Abutment Shoulder 3.7 mm	10
AT42	Abutment Shoulder 4.2 mm	10
EAAx	EAx Analogue line	10
KMP	Manual Adapter for Ratchet Instruments	2
62T	Torx wrench with ratchet attachment	2
CB	Cutter Bar	2
Tray	Instrument Caddy	1

Products available to purchase separately	
Code	Description
CM1	Male Clip - Cylindrical Hole 1 mm Thickness
CM2	Male Clip - Cylindrical Hole 2 mm Thickness
CM3	Male Clip - Elliptical Hole 2 mm Thickness
CF1	Female Clip
BT0	Bar without Clip L 30 mm divisible
BT1	Bar + Male Clip Cylindrical Hole 1 mm thickness
BT2	Bar + Male Clip Cylindrical Hole 2 mm thickness
BT3	Bar + Male Clip Elliptical Hole 2 mm thickness
AT17	Abutment Shoulder 1.7 mm
AT22	Abutment Shoulder 2.2 mm
AT27	Abutment Shoulder 2.7 mm
AT32	Abutment Shoulder 3.2 mm
AT37	Abutment Shoulder 3.7 mm
AT42	Abutment Shoulder 4.2 mm
EAAx	EAx Analogue line
KMP	Manual Adapter for Ratchet Instruments
62T	Torx wrench with ratchet attachment
CB	Cutter Bar
Tray	Instrument Caddy

... Digital  
WINSIX

# WINScan Abutments

# WINScan Abutments

Product Code

Product

**WINScan ABUTMENTS in PEEK**  
Impression or Digital Scan



WS33  
WS33X  
  
WSEA  
WSEAX

**Free Lock Scan Abutment**  
Free Lock Scan Abutment for external hex

**Scan Abutment for EA**  
Scan Abutment for EAx for external hex

**WINBASE ROTATIONAL ABUTMENTS**  
in titanium, for CAD CAM multi-unit screw-retained  
protheses internal hex implants  
VP2 prosthetic screw included



WB3300R  
WB3800R  
WB4500R

● 3.3 mm.  
● 3.8 mm.  
● 4.5 mm.

**WINBASE ROTATIONAL ABUTMENTS**  
in titanium, for CAD CAM multi-unit screw-retained  
protheses external hex implants  
VP2x prosthetic screw included



WB3300Rx  
WB3800Rx  
WB4500Rx

● 3.3 mm.  
● 3.8 mm.  
● 4.5 mm.

Product Code

Product

**WINBASE ANTIROTATIONAL ABUTMENTS**  
in titanium, for CAD CAM cement and screw retained  
prothesis internal hex implants  
VP2 prosthetic screw included



WB3300  
WB3800  
WB4500

● 3.3 mm.  
● 3.8 mm.  
● 4.5 mm.

**WINBASE ROTATIONAL ABUTMENTS**  
in titanium, for CAD CAM multi-unit screw-retained  
prothesis external hex implants  
VP2x prosthetic screw included



WB3300x  
WB3800x  
WB4500x

● 3.3 mm.  
● 3.8 mm.  
● 4.5 mm.

**WINBASE ROTATIONAL ABUTMENTS**  
in titanium, for CAD CAM multi-unit screw-retained  
prothesis conical connection implants  
prosthetic screw VP2c included



WB3802c  
WB3803c  
WB3804c  
WB3805c  
WB3807c

● 3.8 mm.  
● 3.8 mm.  
● 3.8 mm.  
● 3.8 mm.  
● 3.8 mm.



# WINSIX Prosthetic Solutions

<b>WINSIX IMPLANTS</b> Free Lock internal hex connection	Single cement-retai- ned crown	Single screw-retai- ned crown	Permanent cement- retained prosthesis	Permanent multi-unit screw- retained prosthesis	Removable prosthesis bar	Removable Prosthesis	Just on 4/6 Technique
Temporary Abutment	•	•					
Milling Abutment	•		•				
Straight Abutment / MBC / Multipurpose / Aesthetic	•		•				
15° Angled Abutment / MBC / Multipurpose / Aesthetic	•		•				
25° Angled Abutment / MBC / Multipurpose / Aesthetic	•		•				
Overcasting Abutment	•	•	•				
Castable Abutment	•	•	•				
Castable Drillable Abutment	•	•	•				
Ball Attachment						•	
Locator Attachment						•	
WINZIR Abutment	•	•	•				
PRO Abutment	•	•	•				
Flat Shift Line				•	•		
EA LINE - Extreme Abutment				•	•		•
WB Abutment	•	•	•				

<b>TTx IMPLANTS</b> External hex connection	Single cement-retai- ned crown	Single screw-retai- ned crown	Permanent cement- retained prosthesis	Permanent multi-unit screw- retained prosthesis	Removable prosthesis bar	Removable Prosthesis	Just on 4/6 Technique
Temporary Abutment	•	•					
Milling Abutment	•		•				
Straight Abutment	•		•				
15° Angled Abutment	•		•				
25° Angled Abutment	•		•				
Overcasting Abutment	•	•	•				
Castable Abutment	•	•	•				
Locator Attachment						•	
Flat Shift Line				•	•		
EA LINE - Extreme Abutment				•	•		•
CAB - Clip Abutment Bar							•

<b>TTc IMPLANTS</b> Conical connection	Single cement-retai- ned crown	Single screw-retai- ned crown	Permanent cement- retained prosthesis	Permanent multi-unit screw- retained prosthesis	Removable prosthesis bar	Removable Prosthesis	Just on 4/6 Technique
Temporary Abutment	•	•					
Milling Abutment	•		•				
Straight / Multipurpose / Aesthetic Abutments	•		•				
15° Angled Abutment	•		•				
25° Angled Abutment	•		•				
Multifunctional Line				•	•		•
Ball Attachment							•

<b>TTc IMPLANTS</b> Morse Taper connection	Single cement-retai- ned crown	Single screw-retai- ned crown	Permanent cement- retained prosthesis	Permanent multi-unit screw- retained prosthesis	Removable prosthesis bar	Removable Prosthesis	Just on 4/6 Technique
Temporary Abutment	•						
Milling Abutment	•		•				
Straight Abutment	•		•				
15° Angled Abutment	•		•				
25° Angled Abutment	•		•				



# marketing tools for your Practice

 Video for Waiting Room: *It's a Practice of high quality!*

 Practice Plaque

 PATIENT INFORMATION brochure:  
- Implant Solutions  
- Immediate Implants with JUST ON 4

 Implant Card

 Portfolio of Implantology images

 Patient dedicated site  
[www.implantologiawinsix.com](http://www.implantologiawinsix.com)



...Other  
Products

## Product Code

## Product



Product Code	Product
<b>WINSIX FLASH IMPLANTS</b> in titanium, MRS surface	
25011LBC-MRS	Flash Implant LBC Ø 2.5 x 11 mm. MRS surface
25013LBC-MRS	Flash Implant LBC Ø 2.5 x 13 mm. MRS surface
62000F	Mounter for Flash Implants - Manual tool for extracting implants from primary packaging and tightening them
932600	Circular scalpel Ø 2.6 mm.
925000	Screwdriver for Flash Implants with handpiece connection
A25	WINSIX Flash A25 Analogue Ø 2.5 mm.
BOR MICRO	Box complete with O-ring for WINSIX Flash
C25	WINSIX Flash C25 Castable
IM25	WINSIX Flash IM25 Abutment impression coping

## Product Code

## Product



Product Code	Product
<b>Surgical steel OSTEOTOMES</b>	
CW1	<b>Titanium scalpel</b>
CWW I	<b>First step precision marking</b>
CWW II	<b>Second precision marking</b>
CW0 III	<b>Implant Compactor 3.3 mm.</b>
CW0 IV	<b>Implant Compactor 3.8 mm.</b>
CW0 V	<b>Implant Compactor 4.5 mm</b>
CW6	<b>Complete kit</b>



SZ1	<b>Concave Truncated-Cone Shaped Osteotome Ø 2.5 mm. - h 4.0 mm.</b>
SZ2	<b>Convex Truncated-Cone Shaped Osteotome Ø 2.5 mm. - h 5.5 mm.</b>
SZ3	<b>Concave Truncated-Cone Shaped Osteotome Ø 2.5 mm. - h 7.0 mm.</b>
SZ4	<b>Convex Truncated-Cone Shaped Osteotome Ø 2.5 mm. - h 8.5 mm.</b>
SZ5	<b>Convex Truncated-Cone Shaped Osteotome Ø 3.1 mm. - h 10 mm.</b>
SZA	<b>Adapter for Cylindrical Osteotomes</b>
SZ33	<b>Concave Cylindrical Osteotome Ø 3.3 mm. h 15 mm.</b>
SZ38	<b>Concave Cylindrical Osteotome Ø 3.8 mm. h 15 mm.</b>
SZ45	<b>Concave Cylindrical Osteotome Ø 4.5 mm. h 15 mm.</b>
SZ52	<b>Concave Cylindrical Osteotome Ø 5.2 - mm. h 15 mm.</b>
SZ59	<b>Concave Cylindrical Osteotome Ø 5.9-mm. h 15 mm.</b>



## Metals

Type	Composition
Titanium Titanium	Ti grade 4 ASTM F67 Ti grade 5 ASTM F132
Stainless Steel	AISI 303 - AISI 630 - ASTM F899
Platinum Alloy	Au 15%, Pd 60%, Pt 24.9%, Ir 0.1% <i>Melting range 1350° - 1460° C</i>

## Synthetic Resins

Type	Composition
Polymethyl Methacrylate	PMMA <i>Melting range between 150-160° C</i>
Polyoxymethylene	POM <i>Melting range 164-168° C</i>

- Detergent used: SEKUSEPT or similar, dilution 1 measure of detergent per litre of water.
- Decontaminate in solution for 15 min.
- Ultrasonic cleaning for 15 min. At 60°C.
- Dry (very important).

**NB:** After a checking for the presence of any residue, the equipment must be bagged and sterilised at 134°C for 45 minutes.

The equipment put in the autoclave must be rinsed and dried very well since the autoclaving process increases the oxidising action of detergents

LOCATOR is a trademark  
registered property of: 



# Main Bibliography

## WINSIX Implant System

**Evaluation of Ultrashort and Longer Implants with Microrough Surfaces: Results of a 24-to-36-Month Prospective Study**  
L. Malchiodi, R. Caricasulo, A. Cucchi, R. Vinci, E. Agliardi, E. Gherlone - *The International Journal of Oral & Maxillofacial Implant* 2017

**The “Alternating Osteotome Technique”:** a surgical approach for combined ridge expansion and sinus floor elevation. A multicentre prospective study with a three-year follow-up  
L. Malchiodi, A. Cucchi, P. Ghensi, R. Caricasulo, P.F.Nocini Biotechnology & Biotechnological Equipment 09 June 2016

**Riabilitazioni estetiche complesse: case report**  
A. Cardarelli, F. Cardarelli, R. Vinci, E. Gherlone  
*Quintessenza Implantologia* 2016 / numero 2

**Studio clinico sulle riabilitazioni a ridotto numero di impianti solidarizzati a carico immediato**  
A. Cardarelli - *Doctor OS* Settembre 2016 XXVII 07

**Combined microcomputer tomography, biomechanical and histomorphometric analysis of the peri- implant bone: a pilot study in minipig model**  
M. Gramanzini, S. Gargiulo, F. Zarone, R. Megna, A. Apexlla, R. Aversa, M. Salvatore, M. Mancini, R. Sorrentino, A. Brunetti  
*Dental Materials* 32 ( 2016) 794-806

**Sealing ability to Staphylococcus aureus a 4 different implant-abutment connections**  
E. F. Cagidiaco, G. Landini, Lucia Pallecchi, E. Gherlone, Marco Ferrari - *Poster Collegio dei Docenti* 2016

**Riabilitazione a carico immediato con numero ridotto di impianti**  
A. Cardarelli, F. Cardarelli, R. Vinci  
Università Vita-Salute San Raffaele, Milano Corso di Laurea Magistrale in Odontoiatria e Protesi Dentaria - IRCCS Ospedale San Raffaele, Milano Dipartimento di Odontoiatria  
*San Raffaele Dental Clinic - Dr.Os*, marzo 2016 - XXVII 03

**Pianificazione implantoprotesica digitale. Progettazione e realizzazione di riabilitazioni “full arch” a carico immediato**  
M. Manacorda, R. Vinci, F. Bova, E. F. Gherlone  
Università Vita-Salute San Raffaele, Milano Corso di Laurea Magistrale in Odontoiatria e Protesi Dentaria - IRCCS Ospedale San Raffaele, Milano Dipartimento di Odontoiatria  
*San Raffaele Dental Clinic - Dr.Os*, gennaio 2016 - XXVII 01

**Digital impressions for fabrication of definitive “all on four” restorations**  
E. F. Gherlone, F. Ferrini, R. Crespi, G. Gastaldi, P. Cappare', *Implant Dentistry* Vol. 24 Number 1 2015

**Implant prosthetic rehabilitation in controlled HIV- positive patients: a prospective longitudinal study with one - year follow - up**  
E. F. Gherlone, P. Cappare', S. Tecco, E. Polizzi, G. Pantaleo, G. Gastaldi, G. Grusovin  
*Clinical Implant Dentistry and Related Research* Accepted for publication with minor revision ( Manuscript ID CID - 14-459) 2015

**Equine and porcine bone substitutes in maxillary sinus augmentation: a histological and immunohistochemical analysis of VEGF expression**  
S. Tetè, V.L. Zizzari, R. Vinci, S. Zara, U. Di Tore, M. Manica, A. Cataldi, C. Mortellaro, A. Piattelli, E.F. Gherlone  
*The Journal of Craniofacial Surgery* 2014 - May;25(3):835-9. doi: 10.1097/SCS.0000000000000679

**Seven – years retrospective analysis of implant – prosthetic treatment with Winsix[Simbolo] Implant System**  
P. Cappare', F. Bova, D. Defilippi, E.M. Polizzi, E.F. Gherlone  
*Minerva Stomatologica* 2014; 63(4, Suppl.1): 241

**Connection between prosthetic overlplant mesostructures and abutments, precision perShapence using different realization techniques – a preliminary study**  
M. Manacorda, D. Defilippi, A. Ligabue, G. Gastaldi, R. Vinci  
*Minerva Stomatologica* 2014; 63 (4, Suppl. 1): 210

**Treatment options in cases of misplacement implant**  
A. Cardarelli, A. Ligabue, D. Defilippi, R. Vinci  
*Minerva Stomatologica* 2014; 63 (4, Suppl. 1): 251

**Digital Planning and Surgery: the ultimate image based procedure for a successful immediate loading rehabilitation**  
Sh. Ehsani, M. Manacorda, P. Cappare', R. Vinci, G. Gastaldi  
*Minerva Stomatologica* 2014; 63 (4, Suppl. 1): 25

**Rehabilitation of atrophic maxilla with a minimally invasive technique**  
F. Ferrini, M. Manacorda, A. Cardarelli, P. Cappare', R. Vinci  
*Minerva Stomatologica* 2014; 63 (4, Suppl. 1): 236

**Maxillary Sinus Augmentation with Autologous and Heterologous Bone Graft: A Clinical and Radiographic Report of Immediate and Delayed Implant Placement**  
M. Santagata, U. Tozzi, G. Tartaro, V. Santillo, C. Giovanni, E. Lamart, A. Ibro, G. Colella  
*Journal of Maxillofacial and Oral Surgery* 2013 - DOI 10.1077/s12663-013-0569-5

**Riabilitazione dei mascellari edentuli: presentazione di una nuova tecnica software assistita**  
M. Manacorda, R. Vinci, F. Bova, M. Nagni, E. F. Gherlone  
*Doctor Os*, Vol. XXIV - maggio 2013

**Terapia Implantare in Pazienti HIV positivi: studio clinico prospettico**  
G. Gori, P. Cappare', A. Ligabue, A. Cardarelli, D. DeFilippi, M. Nagni, R. Vinci - *Dental School (Dean Prof. Enrico F. Gherlone)*  
Vita Salute University San Raffaele, Milano - Italy  
Sessione Poster - Collegio dei Docenti Roma 16-20 aprile 2013

**An In Vitro Evaluation of Heat Production during Osteotomy preparation for Dental Implants with compressive Osteotomes**  
A. Quaranta, S. Andrea, L. Spazzafumo, M. Piemontese  
*Implant Dentistry* Vol.22 N°2 2013

**Airpolishing: protocollo sperimentale di studio in vitro per la gestione dei pazienti implantari affetti da mucosite e perimplantite**  
P.P. Paganin - *L'igienista moderno* Marzo 2013

**Maxillary sinus augmentation procedures through equine-derived biomaterial or calvaria autologous bone: immunohistochemical evaluation of OPG/RANKL in humans**  
S. Tetè, R. Vinci, V.L. Zizzari, S. Zara, V. La Scala, A. Cataldi, E.F. Gherlone, A. Piattelli  
*European Journal of Histochemistry* Vol.57:e10 2013

**Crestal bone remodeling around platform switched, immediately loaded implants placed in sites of previous failures**  
A. Quaranta, A. Cicconetti, L. Battaglia, M. Piemontese, G. Pompa, I. Vozza - *European Journal of Inflammation* Vol.10 N°2 Accepted July 30 2012

**Use of Piezosurgery During Maxillary Sinus Elevation: Clinical Results of 40 Consecutive Cases**  
M. Cassetta, L. Ricci, G. Iezzi, S. Calasso, A. Piattelli, V. Perrotti  
*The International Journal of Periodontics & Restorative Dentistry* Vol.32 N°6 2012

**Macroscopic and Microscopic evaluation of a new implant design supporting immediately loaded full arch rehabilitation**  
S. Tetè, V. Zizzari, A. De Carlo, B. Sinjari, E. F. Gherlone  
*Annali di Stomatologia* Vol. 3 n°2 2012

**Injectable Magnesium - Enriched Hydroxyapatite Putty in Peri-Implant Defects: A Histomorphometric Analysis in Pigs**  
R. Crespi, P. Cappare', A. Addis, E. F. Gherlone  
*The International Journal of Oral & Maxillofacial Implants* Vol 27, N° 1 2012

**TEETH JUST ON SIX: Riabilitazione Implantoprotesica nei casi di Edentulia completa: un approccio diagnostico-terapeutico integrato**  
M. Manacorda, R. Vinci, M. Pasi, E. F. Gherlone  
*Odontoiatria, Rivista Amici di Brugg n°1*, anno 30°- Maggio 2011

**Long-Term Evaluation of Maxillary Reconstruction by Iliac Crest Bone Graft: A Morphologic and Immunohistochemical Study**  
R. Vinci, R. Vinci, S. Zara, V.L. Zizzari, A. De Carlo, G. Falco, D. Tripodi, A. Cataldi, C. Mortellaro, E.F. Gherlone  
*The Journal of Craniofacial Surgery* Vol. 22, N°5 – Sept. 2011

**Immunohistochemical analysis of Matrix Metalloproteinase -9, Vascular Endothelial Growth Factor, Bone Sialoprotein and i-Nitric Oxide Synthase in Calvaria vs. iliac Crest Bone Grafts**  
S. Tetè, S. Zara, V.L. Zizzari, A. De Carlo, R. Vinci, A. Cataldi, E. F. Gherlone - *Vol. 23, 11 2012*

**Microcomputed and Histologic Evaluation of Calvarial Bone Grafts: A pilot Study in Humans**  
R. Vinci, A. Rebaudi, P. Cappare', E. F. Gherlone  
*The International Journal of Periodontics & Restorative Dentistry* Vol.31, N° 4 2011

**Tecniche Chirurgiche Piezoelettriche**  
R. Vinci, M. Pasi, S. Tetè, E. F. Gherlone  
*Panorama Dental – Sept. 2011*

**Osteosynthesis plates, screws, xenogenic graft and resorbable barriers for preimplant and peri- implant surgery**  
G. M. Rauso, N. Nesi, L. Fragola, M. Santagata, V. Santillo, R. Rauso - *Minerva Stomatologica – 2010, 59:00-00*

**Utilizzo del Mucotome per il rialzo del seno mascellare con approccio crestale: Report preliminare**  
F. Vannini, M. Nardone - *Italian Oral Surgery* 2010; 9 (6):1-14

**Lembo di bolla di Bichat associato al lembo di Rehrmann nella correzione plastica di ampia fistola oro antrale in paziente trattata con bifosfonati: Case Report**  
M. Del Brutto, G. Alfieri, G. Mariani, C. Vianale, R. Mazzanti, G. Sammartino - *Sessione Poster- Roma* 18-19-20 Febbraio 2010

**Difetti di spessore crestali: Utilizzo della tecnica Split Crest**  
B. Marcelli, A. Delli Carpini, E. Coccia, M. Proccaccini  
Sessione Poster-Roma 18-19-20 Febbraio 2010

**“Influence of Novel Nano-Titanium Implant Surface on Human Osteoblast Behavior and Growth”**  
S. Tetè, F. Mastrangelo, R. Quaresima, R. Vinci, G. Sammartino, L. Stuppia, E. F. Gherlone, - *Implant Dentistry*, Vol.19, N° 6 2010

**“Chiave dinamometrica digitale: applicazioni in Implantoprotesi”**  
L. Prosper, M. Cappello, C. A. Cortella, T. D’Alicandro, E. F. Gherlone  
*QI Quintessenza Internazionale*, anno 26 n. 1 Gennaio/Marzo 2010

**“L’impiego della teca cranica a scopo implantoprotesico”**  
E. F. Gherlone, R. Vinci, L. D’Aversa  
*Oral & Implantology* Anno II N°2/ 2009

**“Effect of Implant Angulation, Connection Length, and Impression Material on the Dimensional Accuracy of Implant Impressions: An In Vitro Comparative Study”**  
R. Sorrentino, E. F. Gherlone, G. Calesini, F. Zarone  
*Clinical implant dentistry and related Research*, 2009

**“A Randomized Prospective Multicenter Trial Evaluating the Platform-Switching Technique for Prevention of Postrestorative Crestal Bone Loss”**  
L. Prosper, S. Redaelli, M. Pasi, F. Zarone, G. Radaelli, E. F. Gherlone - *JOMI The international Journal of Oral & Maxillofacial Implants – vol. 24 n. 2, 2009*

**“Isolation of osteogenic progenitors from human amniotic fluid using a single step culture protocol”**  
I. Antonucci, I. Iezzi, E. Morizio, F. Mastrangelo, A. Pantalone, M. Mattioli Belmonte, A. Gigante, V. Salini, G. Calabrese, S. Tetè, G. Palka, L. Stuppia  
*SILENCE a Journal of RNA regulation*, 2009

**“Experimental in vitro study for the implementation of air polishing during treatment of implanted surfaces in patients with mucositis/peri-implantitis”**  
P. Paganin, C. Cortella, E. M. Polizzi, E. F. Gherlone  
*Europerio 6th*, June 4-6 2009 Stoccolma, Svezia

**“In vitro behaviour onto different titanium surface of osteoblast-like cells obtained from human dental pulp”**  
S. Tetè, F. Mastrangelo, V. Zizzari, G. D’Apolito, N. Fiorentino, U. Desiato, M.T. Sberna, R. Quaresima, L. Stuppia, R. Vinci, E. F. Gherlone - *Atti del 7th Annual Meeting of ISSCR Int. Society of Stem Cell Research, Barcellona* July 2009

**“Novel Protocol of osteogenic differentiation from amniotic fluid cells”**  
S. Tetè, F. Mastrangelo, M. Tranasi, V. Zizzari, I. Antonucci, G. D’Apolito, T. Marchese, R. Vinci, L. Stuppia, E. F. Gherlone  
*Atti I Expo di Autunno Università della Lombardia, 1 Dicembre 2007*

**“Interfaccia osso- Implant nei differenti tipi di carico degli impianti dentali”**  
S. Tetè, G. D’Apolito, F. Mastrangelo, R. Vinci, E. F. Gherlone  
*Atti III Expo di Autunno Università della Lombardia, 27-28 Novembre 2009*

**“Valutazione della capacità osteogenica di hafscs ottenute da liquido amniotico”**  
S. Tetè, U. Di Tore, V. Zizzari, L. Stuppia, F. Zarone, E. F. Gherlone  
*Atti III Expo di Autunno Università della Lombardia, 27-28 Novembre 2009*

**“Analisi Morfostrutturale di Pamcs umane sulla Surface a deposizione elettrochimica FCC”**  
S. Tetè, A. De Carlo, E. Nargi, V.F. Zizzari, I. D’Alimonte, R. Ciccarelli, R. Vinci, E. F. Gherlone  
*Atti III Expo di Autunno Università della Lombardia, 27-28 Novembre 2009*

**L’innesto onlay autologo nelle atrofie localizzate**  
B. Marcelli - *QI Quintessenza Internazionale*, anno 25 numero X

**“Evaluation of Effects on bone Tissue of Different Osteotomy Techniques”**  
S. Tetè, R. Vinci, V. Zizzari, L. Cingano, R. Bollero, G. D’Apolito, L. Benzo, A. Dolci, C. Mortellaro, E.F. Gherlone  
*The Journal of Craniofacial Surgery* vol. 20 n. 5 September 2009

**“Osteotomie eseguite con Piezosurgery: analisi biomolecolare degli effetti sul tessuto osseo”**  
S. Tetè, R. Vinci, F. Mastrangelo, U. Desiato, V. Zizzari, N. Fiorentino, L. Renzo, E. F. Gherlone  
*Atti III Expo di Autunno Università della Lombardia, 27-28 Novembre 2009*

**La chirurgia Orale piezoelettrica**  
M. Pasi, R. Vinci, D. Di Stefano, E. F. Gherlone  
*IOS Italian Oral Surgery* 5/2008

**“In vitro evaluation of osteoblast-like cells from different sources”**  
F. Mastrangelo, M. Tranasi, V. Zizzari, D. Farronato, T. Traini, R. F. Grassi, L. Stuppia, S. Tetè  
*Atti del 86th General Session of International Association of Dental Research ( I.A.D.R.) Toronto, 2-5 Giugno 2008*

**“A macro- and nanostructure evaluation of a Novel Dental Implant”**  
S. Tetè, F. Mastrangelo, T. Traini, R. Vinci, G. Sammartino, G. Marezzi, E. F. Gherlone - *Implant Dentistry* vol.17n.3 2008

**“Calvarial versus Iliac Crest for Autologous bone graft material for a Sinus lift procedure: a Histomorphometric Study”**  
R. Crespi, R. Vinci, P. Cappare', E. F. Gherlone, G. E. Romanos  
*JOMI The International Journal of Oral & Maxillofacial Implants*, vol. 22 n. 4 a. 2007

**“A Randomized Prospective Multicenter Trial Evaluating the Platform-Switching Technique for Preventing Postrestorative Crestal Bone Loss”**  
L. Prosper, S. Redaelli, M. Pasi, F. Zarone, G. Radaelli, E. F. Gherlone

**Le tecniche piezoelettriche, recentemente impiegate in chirurgia dei mascellari, consentono osteotomie precise e sicure nel rispetto delle strutture ossee**  
R. Vinci - *Dental Cube* 2 -Marzo 2007

**“Valutazione del comportamento di osteoblasti derivanti da cellule staminali di liquido amniotico su differenti superfici implantari”**  
A. Desiderio, A. D’Incecco, L. Stuppia, I. Antonucci, S. Tetè  
*Atti I Expo di Autunno Università della Lombardia, 1 Dicembre 2007*

**“Evaluation of surface nano-topographic effect on afscells growth”**  
A. Desiderio, F. Mastrangelo, T. Traini, L. Stuppia, G. Sammartino, S. Tetè - *Atti dell’Annual Meeting del I.A.D.R. International Association of Dental Research, Continental Division Thessaloniki* 26-29 Settembre 2007

**“Implant troncoconico con componente protesica sottodimensionata per mantenere l’osso crestale”**  
L. Prosper, S. Redaelli, A. D’Addona, E. F. Gherlone  
*Poster 13° Congresso Nazionale del Collegio dei Docenti di Odontoiatria, Roma* 5-8 Aprile 2006

**“Passivazione della protesi impianto-supportata”**  
L. Prosper, S. Redaelli, T. D’Alicandro, A. D’Addona, E. F. Gherlone  
*Implantologia Orale: rivista di clinica, ricerca, qualità - anno 8, n. 5 2005*

**“Stabilizzazione della protesi mediante ancoraggio con mini impianti in pazienti con edentulia totale inferiore e atrofia del processo alveolare”**  
D. Perazzo, E. Cavalieri, S. Piovani, E. Stellini, E. D’Andrea  
*Poster 12° Congresso Nazionale Collegio dei Docenti di Odontoiatria, Roma - Palazzo dei Congressi* 16 -19 Marzo 2005

**“Follow-up a 4 anni di impianti post-estrattivi a diametro maggiorato”**  
L. Prosper, S. Redaelli, A. D’Addona, E. F. Gherlone  
*Poster n° 486 11° Congresso Nazionale - Collegio dei Docenti di Odontoiatria, 2004*

**“Follow-up a 4 anni di impianti con diametro maggiore inseriti in alveoli di recente estrazione in combinazione con una membrana riassorbibile o con un materiale alloplastico riassorbibile”**  
L. Prosper, E. F. Gherlone, S. Redaelli, M. Quaranta  
*QE Implantologia la Rivista per il Clinico* anno 2- Aprile 2004

**“4 Year follow-up of larger diameter implants placed in fresh extraction sockets using a resorbable membrane or resorbable alloplastic material”**  
L. Prosper, E. F. Gherlone, S. Redaelli, M. Quaranta  
*JOMI The International Journal of Oral & Maxillofacial Implants*, vol. 18 - n. 6 - 2003

**“A system for the diagnosis, placement, and prosthetic restoration of root form implants” (U.S. Patent # 5,769,636)**  
F. Di Sarò  
*Journal of PROSTHODONTICS* , vol. 12 – n. 1 March - 2003

**“Sistema implantare WINSIX : Case-Report di un monoplant post-estrattivo in zona estetica”**  
L. Prosper  
*QO Quintessenza Odontotecnica*, anno 19, Marzo 2002 Verona

**“Osteointegrazione a 5 anni di impianti sabbiati o sabbiati e mordenzati ritenenti protesi parziali fisse”**  
L. Prosper, F. Di Carlo, I. Vozza, M. Quaranta  
*QI Quintessenza Internazionale* I - 2002

**“La passivazione implantare diretta secondo la metodica Free Tense TM” - Winsix Ltd, London, UK**  
S. Redaelli, L. Prosper, E. F. Gherlone, A. D’Addona - 2002

**“Il carico immediato: studio prospettico e comparativo tra due diverse metodiche terapeutiche”**  
S. Redaelli - 2001

**“Clinical trial on osseointegration using sandblasted or sandblasted and acidetched implants”**  
S. Redaelli, L. Prosper, F. Di Carlo, A. Daddona, M. Quaranta - 2001

**“Implant-retained mandibolar overdentures with immediate loading: clinical and histological study”**  
G. Pompa, M. Cassetta, I. Vozza, R. Scaringi, M. Quaranta  
*JDR, Journal of Dental research* vol. 80, March 2001

**“ Healing period of titanium implants with sandblasted and acid-etched surface”**  
L. Prosper, F. Di Carlo, S. Redaelli, R. Scaringi, M. Quaranta  
*3° World Congress of Osteointegration* 2001

**“Tempo di guarigione ossea per impianti di titanio con Surface sabbiata o Surface sabbiata e mordenzata”**  
L. Prosper, F. Di Carlo, S. Redaelli, G. Radaelli, R. Scaringi, M. Quaranta - *La rivista Internazionale di Odontoiatria Protesica*, vol.13 n. 1 - 2000

**“Using of Gore Resolut Membrane in 223 post-extraction WINSIX Implants”**  
L. Prosper, S. Redaelli, G. Pompa, A. Palattella, M. Quaranta  
*Atti 4TH Joint meeting* 24-27 August 2000, Varsavia - Polonia

**“Gore Resolut membrane and implantation below the crest in postextraction surgery”**  
F. Di Carlo, L. Prosper, S. Redaelli, I. Guadagno, M. Quaranta – anno 2000

**“Analisi in vivo su tre superfici implantari: valutazione istomorfometrica”**  
S. Redaelli, L. Prosper, F. Di Carlo, A. Scaranò, I. Vozza – anno 2000







**BIOSAFIN srl**  
info@biosafin.com  
biosafin.com

OFFICES IN ITALY

**MILAN:** Via Cagliari 32/44  
20060 - Trezzano Rosa - Zona Industriale (MI)  
Tel. +39 02 90968692  
Fax +39 02 90968541

**ANCONA:** Via Tiraboschi, 36/G  
60131 - Ancona (AN)  
Tel. +39 071 2071897  
fax +39 071 203261

# BIOSAFIN



**BIOSAFIN is a Certified Company:**

**UNI EN ISO 9001:2015**, which certifies the entire 360° work process, guaranteeing compliance with the quality standards considered to be optimum for the protection of the Product User - the Professional - and the end user - the Patient.

**UNI CEI EN ISO 13485: 2016** specifically relevant to the Quality of Medical Devices.

**WINSIX**<sup>®</sup>  
PERFORMING IMPLANT SYSTEM 

PRODUCTS CERTIFICATIONS  
**CE** **FDA**