

HaemoCer™ PLUS - versatile use:

- Treatment of diffuse and localized bleeding
- For bleeding around suture lines
- Substitution of fibrin glue for haemostasis in vascular surgery
- Substitution of bone wax in bone bleeding
- Use with compatible filter cell-savers

Fields of application:

- Cardiac, Vascular and Thoracic Surgery
- General Surgery
- Transplant Surgery
- Traumatology
- Tumour Surgery
- Oral and maxillofacial Surgery
- ENT Surgery
- Urology / Gynaecology
- Orthopaedics
- Burn surgery
- Neurosurgery

HaemoCer™ PLUS

Haemostatic powder including prophylactic functions preventing adhesions¹⁾

Product Overview HaemoCer™ PLUS

HFP201	1g	5 pcs
HFP202	2g	5 pcs
HFP203	3g	5 pcs
HFP205	5g	5 pcs
HFZ101	440 mm	5 pcs
HFZ201	180 mm	5 pcs



BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0
 Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de



HDD201_10/2020-09

Natural - Safe - Effective

BioCer

HaemoCer™ PLUS

HaemoCer™ PLUS Absorbable Polysaccharide Haemostat (APH) is a proprietary patent pending technology based on BioCers *Polysaccharide Ultra-hydrophilic Resorbable Engineering (PURE)* process. HaemoCer™ PLUS is a resorbable plant based haemostat containing no

animal nor human components. Haemostasis occurs by rapidly accelerating the normal physiologic clotting cascade without participation of chemical or pharmaceutical ingredients.

Natural

- Plant based
- Contains no animal nor human components
- Completely absorbable

Safe

- Rapid resorption minimizes risks of foreign body reactions and granuloma formation
- Reduction of homologous transfusions and blood donations

Effective

- Double enhanced absorption capacity, compact packaging reducing storage volume
- Rapid off-the-shelf deployment, no special storage, defrosting or premixing

Ease of use, DAPI application technique

DRY

Remove excess pooled blood by suction or gauze



APPLY

Apply HaemoCer™ PLUS immediately and extensively to the wound site



PRESS

Apply gentle compression over the wound using dry gauze



IRRIGATE

Irrigate thoroughly prior to gently removing the gauze to ensure the gelled matrix remains at the wound site



Nature refined

HaemoCer™ PLUS – Absorbable Polysaccharide Haemostat (APH) engineering incorporates a sophisticated, plant-based polymer crosslinking that creates ultra-hydrophilic, biocompatible particles.

Once in contact with blood, HaemoCer™ PLUS enhances the natural clotting cascade by rapidly dehydrating the blood and accelerating the concentration of platelets, red blood cells and coagulation proteins at the bleeding site.

The second mode of action of HaemoCer™ PLUS is the formation of a robust gelled matrix that adheres to the bleeding site and forms a mechanical barrier to prevent further bleeding.

1) Studies with both HaemoCer™ PLUS Haemostatic Powder and comparable competitor

products have shown that the gel acts as a mechanical barrier and helps reduce postoperative adhesions.

HaemoCer™ PLUS contains no animal nor human components and is rapidly resorbed by amylase within a few days without leaving any residues behind, depending on the amount of used material and the place of application.

Plant Based

