Implanted **electronic devices** (for example a cardiac pacemaker) - specialist medical opinion should be obtained before the patient is subjected to stimulation. Recently diagnosed or suspected DVT: specialist medical opinion should be obtained before the patient is subjected to stimulation.

Read the instructions carefully before use. Only apply the device in accordance with instructions. Improper use of this device can be dangerous.

The geko™ device may be operated either by the patient or a third party. Patients should be trained in the positioning and operation of the device.

Keep out of the reach of children and pets. No modification of this equipment is allowed. Do not use whilst driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user or others at undue risk of injury.

Stimulation should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck, (especially the carotid sinus), or from electrodes placed on the chest and the upper back or crossing over the heart.

Application of electrodes near the thorax may increase the risk of cardiac fibrillation.

Do not apply over sore, infected or inflamed areas, broken skin or skin eruptions, ecz, phlebitis, thrombophlebitis, varicose veins etc. Do not apply over, or in proximity to, cancerous lesions. Do not use in proximity (i.e. within 1m) of short wave/microwave equipment, as this may affect the device.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the geko. Otherwise, degradation of the performance of this equipment could result.

Do not use if packaging is open or appears damaged. Do not use the device that appears to be damaged. Do not use the device in close proximity to water sources, such as fires or radiant heaters, as excess heat may be disposed of with care. Batteries must be disposed of in accordance with any local legislation. Some hospital and clinics will have specific requirements for disposal of used medical devices. If used at home the device may be disposed of in your general waste if regulations permit.

For ease of disposal the electrodes may be removed with scissors if necessary. The battery can be removed by breaking open the housing and prising out. See www.gekodevices.com for guidance.

Do not incinerate the device.

The hydrogel contains acrylic polymers. Do not use on patients known to be allergic to acrylic acid. Medical advice should be sought when the patient is pregnant, has diagnosed heart conditions or epilepsy, or following surgery where muscle contractions may disrupt the healing process.

The device should be removed before the patient undergoes MRI as it contains ferromagnetic components.

The device should be switched off during ECG monitoring using leg electrodes as it may interfere with ECG leg electrode signals.

The device is not intended for use in an oxygen rich environment.

Do not place the device in the mouth. The device is a choke and potential allergic hazard. If the device or

any component is swallowed IMMEDIATE medical assistance. Serious harm could be caused if the battery is swallowed.

Check the device is functioning correctly every 2 hours.

Use the device only as instructed. Visit www.gekodevices.com for further advice and usage tips.

The safety of the geko has not been tested on children.

Do not switch on unless correctly attached to the patient.

The device must be kept clean and dry. Do not bath or shower whilst wearing the device.

Do not apply a device that has been worn by another individual as this will carry a risk of infection.

Specialist medical opinion should be obtained before devices are used for more than 28 consecutive days.

Remove the device carefully, lifting the long-end first and applying a constant gentle force.

Excessive force may damage the device.

EMC Declaration

Use of this equipment adjacent to other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

The geko™ device is certified to EN 60601-1-2:2015 regarding Electromagnetic Compatibility. Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided at: **http://www.gekodevices.com/**

The geko may switch off if exposed to high levels of electromagnetic disturbance.

Disposal

The device does not contain any toxic or environmentally hazardous materials. After use the device may be potentially contaminated or infected because it has been in contact with skin for several hours, and so need to be disposed of with care. Batteries must be disposed of in accordance with any local legislation. Some hospital and clinics will have specific requirements for disposal of used medical devices. If used at home the device may be disposed of in your general waste if regulations permit.

For ease of disposal the electrodes may be removed with scissors if necessary. The battery can be removed by breaking open the housing and prising out. See www.gekodevices.com for guidance.

Do not incinerate the device.

FirstKind Limited*

Hawk House Peregrine Business Park Gomm Road High Wycombe HP13 7DL United Kingdom

T: +44 (0)1494 572040 www.gekodevices.com

If assistance is required or to report difficulties, unexpected operation or events, contact the manufacturer.

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