

**CURRENTBODY SKIN™**






# CURRENTBODY

THE BEAUTY TECH EXPERTS

We are the Beauty Tech Experts. Over the past 14 years, we have made it our mission to source breakthrough, clinically proven technologies and bring them to users in the comfort of their own homes.

Home-use beauty tech can give you skincare results beyond the effects of any topical serums and creams alone. **Beauty technology has the power to go deeper than the surface of the skin.**





# CURRENTBODY SKIN™

Developed through our expertise in beauty technology, CurrentBody Skin is a range of high-performance beauty devices for home use, including our world-leading **LED Light Therapy Face Mask**, which has been used by over 500,000 people across 80 countries.

Each of our devices pairs state-of-the-art, clinically-proven technologies, such as **LED Light Therapy and Radio Frequency**, with unique software to give professional-level results at home.

# GLOBAL PRESENCE






# GLOBAL RETAIL PARTNERS



John Bell & Croyden | London



Harrods | London



Selfridges | London



# GLOBAL RETAIL PARTNERS



Printemps | Paris



La Samaritaine | Paris



El Corte Inglés | Spain



# GLOBAL RETAIL PARTNERS



Joyce Beauty | Hong Kong and China




Stockmann | Finland



Kadewe | Germany



# LOVED BY GLOBAL CELEBRITIES



58 MILLION VIEWS ON  
**EMILY  
IN PARIS**




LAUREN RIDLOFF



CILLIAN MURPHY




WAYNE GOSS



OLIVIA DEJONGE



SUKI WATERHOUSE



HALLE BERRY



KIM KARDASHIAN



LISA SNOWDON



NAOMI WATTS



RENÉE ZELLWEGER



KALEY CUOCO



TRINNY WOODALL

GLAMOUR GRAZIA COSMOPOLITAN ELLE THE TIMES



# CURRENTBODY SKIN™



## LED 4-IN-1 ZONE FACIAL MAPPING MASK

109% more powerful than other multi-treatment LED masks



## LED LIGHT THERAPY MASK

The original anti-ageing, flexible LED face mask.



## LED LIGHT THERAPY FACE MASK: SERIES 2 *Coming Soon*

Upgraded light coverage provides our famous instant glow, for a healthier, brighter and smoother complexion.



## LED HAIR GROWTH HELEMT

For advanced hair growth at home.



## ANTI-BLEMISH MASK

Treats mild to moderate acne through blue and red light therapy.



## LED NECK & DEC PERFECTOR

The secret to Kourtney Kardashian's Oscar glow.



CLINICALLY PROVEN



APPROVED

Premium but accessible | RRP between: \$79- \$859



# CURRENTBODY SKIN™



## LED HAND PERFECTOR

Take care of your anti-ageing concerns with a daily 10-minute sessions.



## LASER HAIR REMOVAL

The only at home laser hair removal treatment for permanent, full body results outside of a clinic



## ANTI-BLEMISH LED PEN

Experience the ultimate solution for on-the-go skincare with our portable, fast-acting blemish treatment.



## LED EYE PERFECTOR

LED anti-ageing eye mask for younger-looking eyes.



## LED LIP PERFECTOR

Youthful-looking lips in 3 minutes a day.



## TEETH WHITENING KIT

Clinically-proven LED solution for whiter teeth and healthier gums



CLINICALLY  
PROVEN



APPROVED



# CURRENTBODY SKIN™



## RF RADIO FREQUENCY SKIN TIGHTENING DEVICE

Tightens loose, sagging skin - 89% saw an improvement in skin tightness in 8 weeks



## DERMALUX FLEX MD LED LIGHT THERAPY DEVICE

Award-winning light therapy device trusted by skincare professionals worldwide

# CURRENTBODY WELLNESS™



## INFRARED SAUNA BLANKET

Detox the body and relax the mind with power of infrared heat therapy.



CLINICALLY  
PROVEN



APPROVED



# CURRENTBODY SKIN™



## HYDROGEL FACE MASK

Improve the performance of your LED devices and instantly hydrate your skin



## HYDROGEL DÉCOLLETAGE MASK

Skincare mask to be used under your CurrentBody Skin LED Neck & Décolletage Mask or on its own



## HYDROGEL EYE MASK

Skin is left feeling hydrated, plumped and nourished after one use



## HYDROGEL LIP MASK

Enhance the performance of your LED device to plump and hydrate lips



## GREEN TEA SERUM

Green tea and vitamin C serum to enhance LED mask benefits



## RADIO FREQUENCY CONDUCTIVE GEL

Conducts and enhances radio frequency energy output of the CurrentBody Skin RF Device





# CURRENTBODY SKIN™

## LED LIGHT THERAPY FACE MASK

### THE BIGGEST BREAKTHROUGH IN BEAUTY

Reduce wrinkles by 24% in just 4 weeks with red & near-infrared light.

97% said skin looked **brighter or fresher** after using the LED Face Mask.

Our Pillow Technology covers all areas of the face for **maximum results**.

Quick, 10-minute treatment that can be integrated into your daily skincare routine.







## MARKET POTENTIAL

Significant margin opportunity.

-


UK average monthly searches for **'led face mask' +83% YoY.**

-

The LED Mask was seen by 58 million households alone on Emily in Paris.

-

International fulfillment center: EU, UK, HK, US and AU.



## BEAUTY EXPERTISE

In-house product training with device specialists through 20 years of industry experience.

-

Expert point of sale marketing materials.

-

In-house product development team.



## CELEBRITY & EXPERT APPROVED


Internationally endorsed by celebrities, beauty experts and editors.

-

Features and awards from The Times, Vogue, and Harper's Bazaar and more.

-

Global A-Lister exposure has propelled our products to a global audience.



## IN HOUSE GLOBAL MARKETING

International PR support - delivered over 30 billion reach in 2023.

-

Ongoing global marketing activity.

-

11 global PR agencies.