



**simpass**

# END018

.018" OTW PTA BALLOON DILATATION CATHETER

## **LOW PROFILE / STRONG CATHETER SHAFT**

OTW catheter designed specially for below the knee

## **LONG BALLOON**

Balloon lengths up to 220 mm enables single dilatation for long lesions

## **CATHETER COATING**

Hydrophilic coated OTW shaft for improved crossability

## **SUPERIOR CROSSING PERFORMANCE**

Balloon folding geometry minimizes profile & allows smooth rewrap

Product Code (80 / 130 / 150 cm)	Balloon Diameter (mm)	Balloon Length (mm)	Nominal / RBP (atm)	Recommended Introducer Sheath (F)
SND20020-xx	2.0	20	8 / 16	4
SND20030-xx	2.0	30	8 / 16	4
SND20040-xx	2.0	40	8 / 16	4
SND20060-xx	2.0	60	8 / 16	4
SND200100-xx	2.0	100	8 / 16	4
SND200150-xx	2.0	150	8 / 16	4
SND200220-xx	2.0	220	8 / 16	4
SND25020-xx	2.5	20	8 / 16	4
SND25030-xx	2.5	30	8 / 16	4
SND25040-xx	2.5	40	8 / 16	4
SND25060-xx	2.5	60	8 / 16	4
SND250100-xx	2.5	100	8 / 16	4
SND250150-xx	2.5	150	8 / 16	4
SND250220-xx	2.5	220	8 / 16	4
SND30020-xx	3.0	20	8 / 16	4
SND30030-xx	3.0	30	8 / 16	4
SND30040-xx	3.0	40	8 / 16	4
SND30060-xx	3.0	60	8 / 16	4
SND300100-xx	3.0	100	8 / 16	4
SND300150-xx	3.0	150	8 / 16	4
SND300220-xx	3.0	220	8 / 16	4
SND35020-xx	3.5	20	8 / 16	4
SND35030-xx	3.5	30	8 / 16	4
SND35040-xx	3.5	40	8 / 16	4
SND35060-xx	3.5	60	8 / 16	4
SND350100-xx	3.5	100	8 / 16	4
SND350150-xx	3.5	150	8 / 16	4
SND350220-xx	3.5	220	8 / 16	4
SND40020-xx	4.0	20	8 / 16	4
SND40030-xx	4.0	30	8 / 16	4
SND40040-xx	4.0	40	8 / 16	4
SND40060-xx	4.0	60	8 / 16	4
SND400100-xx	4.0	100	8 / 16	4
SND400150-xx	4.0	150	8 / 16	4
SND400220-xx	4.0	220	8 / 16	4
SND50020-xx	5.0	20	8 / 16	5
SND50030-xx	5.0	30	8 / 16	5
SND50040-xx	5.0	40	8 / 16	5
SND50060-xx	5.0	60	8 / 16	5
SND500100-xx	5.0	100	8 / 16	5
SND500150-xx	5.0	150	8 / 16	5

Product Code (80 / 130 / 150 cm)	Balloon Diameter (mm)	Balloon Length (mm)	Nominal / RBP (atm)	Recommended Introducer Sheath (F)
SND500220-xx	5.0	220	8 / 16	5
SND60020-xx	6.0	20	8 / 16	5
SND60030-xx	6.0	30	8 / 16	5
SND60040-xx	6.0	40	8 / 16	5
SND60060-xx	6.0	60	8 / 16	5
SND600100-xx	6.0	100	8 / 16	5
SND600150-xx	6.0	150	8 / 16	5
SND600220-xx	6.0	220	8 / 16	5
SND70020-xx	7.0	20	8 / 16	5
SND70030-xx	7.0	30	8 / 16	5
SND70040-xx	7.0	40	8 / 16	5
SND70060-xx	7.0	60	8 / 16	5
SND700100-xx	7.0	100	8 / 16	5
SND700150-xx	7.0	150	8 / 16	5
SND700220-xx	7.0	220	8 / 16	5
SND80020-xx	8.0	20	8 / 16	6
SND80030-xx	8.0	30	8 / 16	6
SND80040-xx	8.0	40	8 / 16	6
SND80060-xx	8.0	60	8 / 16	6
SND800100-xx	8.0	100	8 / 16	6
SND800150-xx	8.0	150	8 / 16	7
SND800220-xx	8.0	220	8 / 16	7

Please consult sales representative for different size requests.

<b>Nominal Pressure</b>	8 atm
<b>Rated Burst Pressure</b>	16 atm
<b>Markers</b>	Swaged 1 mm Platinum Iridium Markers (Distal & Proximal)
<b>Guidewire Compatibility</b>	0.018"
<b>Tip Entry Profile</b>	0.022"
<b>Sheath Compatibility</b>	2 mm - 4 mm : 4 F 5 mm - 7 mm : 5 F 8 mm : 6 F 8 mm : 7 F (for 120 mm and longer balloons)
<b>Balloon Fold</b>	2 mm - 6 mm : 3 Fold 7 mm - 8 mm : 6 Fold
<b>Catheter Coating</b>	Hydrophilic Coating

REV04-140222-B