

CovaView™ IV

Transparent IV Dressing



✓ **EXCELLENT MOISTURE VAPOUR TRANSMISSION RATE (MVTR)**

Moisture and oxygen vapour exchange allows for effective adhesion and prevents skin maceration.

✓ **PROTECTS I.V. SITE**

Provides a waterproof, sterile barrier to external contaminants including liquids, bacteria and viruses.

✓ **I.V. SITE VISIBILITY**

Transparent dressing allows for continuous visibility to the I.V. insertion site.

✓ **EASY TO APPLY**

Window frame delivery makes accurate placement easy.

WHAT IS IT?

CovaView™ IV is a thin, polyurethane film dressing with release paper (frame-style). CovaView™ IV is waterproof, and acts as a barrier against viruses, bacteria, and outside contaminants. Its permeable semi-occlusive PU film is breathable, allowing good oxygen and moisture vapour exchange. CovaView™ IV is designed to be applied over I.V. sites, clean skin, or wounds.

WHEN DO I USE IT?

CovaView™ IV is indicated to cover and protect I.V. sites, surgical incisions, primary dressing on pressure ulcers (stages I and II) with minimal drainage, partial thickness wounds, and skin tears or abrasions.

Please refer to package insert for full instructions.

AVAILABILITY:

COVAVIEW™ IV			
Transparent IV Dressing			
Reference Number		Size (cm)	Dressings per Carton Cartons per Case
6000607 Cover Dressing		6cm x 7cm	50/24
6001012 Cover Dressing		10cm x 12cm	50/24
6030785		7cm x 8.5cm	50/24
6038510		8.5cm x 10.5cm	50/24



CovaView™ IV
Transparent IV Dressing

1-877-711-6055 | info@covalon.com

