



Cadwell C2 Cart

Schedule a Cascade PRO demo to experience our time-tested IONM solution.

- 16 or 32-channels (buy now or upgrade later)
- Dynamic montaging
- All-modality monitoring
- Remote monitoring

Cascade PRO supports:

- Auditory and visual stimulation
- Electrosurgery and noise detection
- Cart-based and portable configurations
- ES-IX 2: High and low current stimulation
- TCS-4/TCS-1000: Transcranial stimulation
- Surge protection




Cascade® Surgical Studio is Cadwell's dynamic solution for all Cascade IONM systems.


Surgical Studio helps you:

Streamline your workflow

- Switch users during a case without interrupting data acquisition
- Maximize trace windows on the fly
- Share and manage case data and resources across your organization with CadLink™ information management




View live and saved data side-by-side




Set up color-changing cursor alerts


Capture critical information



Auto-save screenshots with every Event Log entry



Save only important EMG data



Choose screenshots to include in reports

Strengthen documentation

- Standardize demographic data collection
- Enter comments just by typing
- Customize Word®-based report templates
- Choose whether to save the Chat log
- Automatically document remote reviewer connections