

rejoin 

Hangzhou Rejoin Mastin Medical Device Co., Ltd

Headquarter Office: 22 Xinyan Rd, Yuhang Economic Development Zone, Hangzhou, Zhejiang, 311100, P.R. of China

Beijing Branch: Room 1605, Jinghui building, No. 118 Jianguo Road, Chaoyang District, Beijing, P.R. of China

Shanghai Branch: 309A, L1-B, Hongqiao World Center, 1588 Zhuguang Rd, Qingpu District, Shanghai, P.R. of China

contact@rejoin.com [@rejoinmedical](https://www.rejoinmedical.com) www.rejoin-medical.com

Internal information is for professional reference only

rejoin 



SHOULDER PRODUCTS

www.rejoin-medical.com



REJOIN Medical



@rejoinmedical



HealFix®



HealFix®

- High-Low double thread design enhances the pull-out strength
- Two independent eyelets provide optimal fixation of both sutures & eliminate suture twining
- Sharpened tip facilitates easy insertion

HealFix Ti Suture Anchor, 2.0mm, w/two #2-0 sutures W&WB	800500600
HealFix Ti Suture Anchor, 2.0mm, w/two #2-0 sutures W&WB, w/needles ..	800500800
HealFix Ti Suture Anchor, 2.8mm, w/#2 suture WB	800501400
HealFix Ti Suture Anchor, 2.8mm, w/two #2-0 sutures W&WB	800501200
HealFix Ti Suture Anchor, 2.8mm, w/two #2-0 sutures W&WB, w/needles ..	800501000
HealFix Ti Suture Anchor, 3.5mm, w/two #2 sutures W&WB	800503500
HealFix Ti Suture Anchor, 3.5mm, w/two #2 sutures W&WB, w/needles	800503300
HealFix Ti Suture Anchor, 5.0mm, w/two #2 sutures W&WB	800503900
HealFix Ti Suture Anchor, 5.0mm, w/two #2 sutures W&WB, w/needles	800503700
HealFix Ti Suture Anchor, 6.5mm, w/two #2 sutures W&WB	800504300
HealFix Ti Suture Anchor, 6.5mm, w/two #2 sutures W&WB, w/needles	800504100

1.6mm Drill for 2.0mm HealFix	801310000
Drill Guide for 2.0mm HealFix	801310100

1.8mm Drill for 2.8mm HealFix	801310200
-------------------------------------	-----------

Drill Guide for 2.8mm/3.5mm HealFix	801310300
---	-----------

2.5mm Drill for 3.5mm/5.0mm HealFix	801310400
---	-----------

3.5mm Drill for 6.5mm HealFix	801310600
-------------------------------------	-----------

TITANIUM SUTURE ANCHOR

HealFix[®] Plus

- Fully helix transitional thread design increases contact areas with cortical bone for high pull-out strength
- Built-in suture hang-up hole eliminates Cheese-Wire effect between sutures and bone after screw implantation
- Self-punching makes implantation easier

HealFix Plus Ti Suture Anchor, 4.5mm, w/two #2 sutures W&WB	800504500
HealFix Plus Ti Suture Anchor, 5.5mm, w/two #2 sutures W&WB	800504500
HealFix Plus Ti Suture Anchor, 6.5mm, w/two #2 sutures W&WB	800504500

3.5mm Drill for 4.5mm HealFix Plus 801310800

4.5mm Drill for 5.5mm HealFix Plus 801310900

5.5mm Drill for 6.5mm HealFix Plus 801311000



HealFix[®] Plus



PEEK SUTURE ANCHOR

FusionLoc[®] Anchor, PEEK

- Open & vented design facilitates bone ingrowth and accelerates healing
- Extended, fully-threaded anchor design provides significantly more threaded engagement, delivering greater pullout strength in poor-quality osteoporotic bone
- Inserter reaches the tip of the anchor, spreads the stress & achieves better torsion force

FusionLoc PEEK Suture Anchor, 4.5mm, w/two #2 sutures WB&W	800601700
FusionLoc PEEK Suture Anchor, 4.5mm, w/#2 suture WB & suture tape W ..	800603300
FusionLoc PEEK Suture Anchor, 5.5mm, w/two #2 sutures WB&W	800605800
FusionLoc PEEK Suture Anchor, 5.5mm, w/#2 suture WB & suture tape W ..	800607300

Threaded Punch for 4.5mm FusionLoc 801311200

Threaded Punch for 5.5mm FusionLoc 801311300



FusionLoc[®] Anchor, PEEK





GripLoc® Anchor, PEEK



GripLoc® Knotless Anchor, PEEK

- Open & vented design facilitates bone ingrowth and accelerates healing
- Extended, fully-threaded anchor design provides significantly more threaded engagement, delivering greater pullout strength in poor-quality osteoporotic bone
- Inserter reaches the tip of the anchor, spreads the stress & achieves better torsion force

GripLoc PEEK Suture Anchor, 3.5mm, w/#2-0 suture W	800610900
GripLoc PEEK Suture Anchor, 4.75mm, w/#2 suture WB	800601700
GripLoc PEEK Suture Anchor, 4.75mm, w/#2 suture W & suture tape WB	800603300
GripLoc PEEK Suture Anchor, 5.5mm, w/#2 suture WB	800605800
GripLoc PEEK Suture Anchor, 5.5mm, w/#2 suture W & suture tape WB	800607300

Punch for 3.5mm GripLoc	801313300
Drill for 3.5mm GripLoc	801313100
Drill Guide for 3.5mm GripLoc	801313200

Punch for 4.75mm/5.5mm GripLoc & 4.5/5.5 OutPrint	801313000
Threaded Punch for 4.75mm GripLoc	801917100
Threaded Punch for 5.5mm GripLoc	801917000



RealLoc® Anchor, PEEK

RealLoc® Anchor, PEEK

RealLoc PEEK Suture Anchor, 2.5mm×8mm	800614100
RealLoc PEEK Suture Anchor, 2.9mm×15.5mm	800614200
RealLoc PEEK Suture Anchor, 4.5mm×24mm	800614400
Punch for 2.5mm RealLoc	RJ8108003
Punch for 4.5mm RealLoc	RJ8108004

1.8mm drill bit for 2.5mm RealLoc	RJ8106001
2.0mm drill bit for 2.5mm RealLoc	RJ8106002
Drill guide for 2.5mm RealLoc	RJ8106003

Drill for 2.9mm RealLoc	RJ8106001
Drill for 2.9mm RealLoc (hard bone)	RJ8106006
Spade Tip Drill for 2.9 mm RealLoc	RJ8106009
Drill Guide, Trocar and Blunt Tip Obturator for 2.9 mm RealLoc	RJ8106008
Drill guide, Trocar Tip Obturator for 2.9 mm RealLoc	RJ8106010
Offset Guide for 2.9 mm RealLoc	RJ8106007



SolidPrint® Anchor, PEEK

SolidPrint® Anchor, PEEK

- Preloaded WilSuture™ Ultra
- Radiolucent and high strength PEEK material
- Ergonomic handle design provides easy & quick release of anchor

SolidPrint PEEK Suture Anchor, 2.3mm, w/#2 suture WB	800610800
Drill Guide for 2.3mm SolidPrint	801311600
Drill for 2.3mm SolidPrint	801311500
Obturator for 2.3mm SolidPrint	801311700

ALL-SUTURE ANCHOR

JoinFix® All-Suture Anchor

- For surgeons that prefer an “all-suture” anchor for the medial row of a double- row, bridging, rotator cuff repair
- 1.5mm WilSutureTape offers: Increased resistance to tissue pull-through. Stronger knotted and knotless fixations. Tighter, smaller knot stacks.
- Deliver secure fixation with up to 80% less bone removal
- UHMWPE is an ideal combination of strength, biocompatibility, revisability and inherent radiolucent characteristics
- Ergonomic handle designed for surgeon comfort and ease of use

JoinFix All-suture Anchor, 1.4mm, w/#2 suture WB	801101300
JoinFix All-suture Anchor, 2.4mm, w/two #2 sutures WB&W	801101800
JoinFix All-suture Anchor, 2.4mm, w/#2 suture W & suture tape WB	801102300
JoinFix All-suture Anchor, 2.4mm, w/#2 suture W & suture tape W	801102400
Drill Guide for 1.4mm JoinFix	801330300
1.6 mm Drill for 1.4mm JoinFix	801913400
Obturator for 1.4mm JoinFix	801330600
Drill Guide for 2.4mm JoinFix	801330500
2.6mm Drill for 2.4mm JoinFix	801903500
Obturator for 2.4mm JoinFix	801330800

JoinFix® Plus

- Radially expanding offers highest fixation strength, ultra-low displacement
- Providing the benefits of a small, soft anchor with the characteristics of traditional anchor designs
- UHMWPE is an ideal combination of strength, biocompatibility, revisability and inherent radiolucent characteristics
- Ergonomic handle designed for surgeon comfort and ease of use

JoinFix Plus All-suture Anchor, 1.7mm, w/#2 suture WB	801100100
JoinFix Plus All-suture Anchor, 2.7mm, w/two #2 sutures WB&W	801100500
JoinFix Plus All-suture Anchor, 2.7mm, w/#2 suture WB & suture tape W ...	801100900
Drill for 1.7mm JoinFix Plus	801330102
Obturator for 1.7mm JoinFix Plus	801330106
Drill Guide for 1.7mm JoinFix Plus	801330108
Kit for 1.7mm JoinFix Plus	RJ8112159
Drill for 2.7mm JoinFix Plus	801330202
Obturator for 2.7mm JoinFix Plus	801330206
Drill Guide for 2.7mm JoinFix Plus	801330208
Kit for 2.7mm JoinFix Plus	RJ8112161



JoinFix®



JoinFix® Plus





ETButton® AC



ETButton® AC

- Four-leaf clover design eases loading 2 Wilsuture® Tapes, leading to stronger flexible fixation
- Large size precontoured button maximizes contact area, reducing the risk of bone collapses
- Laser mark for accurate button loading

ETButton AC 800900600

AC Sleeve Handle 801902000

AC Drill Sleeve, 3mm 801902100

AC Target Guide, left 801902500

AC Target Guide, right 801902600

AC Guide Handle 801902700

AC Cannulated Drill, 2.4mm, Outer 801902200

AC Cannulated Drill, 2.4mm, Inner 801902300

AC Cannulated Drill, 3.0mm 801902400

ACROMIOCLAVICULAR JOINT REPAIR

WilSuture[®] Tape

- Multiple WilSuture[®] Tapes feature outstanding tensile strength for AC joint reduction
- Silicon coating makes knot-tying easy and secures knots

WilSuture Tape Pure, UHMWPE, 1.5mm, 38", white-blue	801243300
WilSuture Tape Pure, UHMWPE, 1.5mm, 38", white	801243400
WilSuture Tape, UHMWPE, 2mm, 7", white-blue	801242700



WilSuture[®] Tape

SHOULDER REPAIR TOOLS

Shoulder Set Box	801302600
Tissue Grasper	801300800
Suture Cutter	801300900
Suture Grasper	801301000
Hook Probe	801301400
Extra Long Switching Stick, 4mm	801301500





Suture Lasso Handle	801302000
Suture Lasso 90° Straight	801302100
Suture Lasso 45° Left	801302200
Suture Lasso 45° Right	801302300
Suture Lasso 90° Left	801302700
Suture Lasso 90° Right	801302800
Suture Lasso Wire Loop	801905900

Shoulder Tissue Elevator, 30°	801301600
-------------------------------------	-----------

Glenoid Rasp	801301700
--------------------	-----------

Knot Pusher, Closed End	801301800
-------------------------------	-----------

Crochet Hook	801301900
--------------------	-----------

EasyPass MultiFire Suture Passer	801302400
--	-----------

EasyPass Needle, MultiFire	801302500
----------------------------------	-----------

EasyPass Suture Passer, Frog 45° Right	801303100
EasyPass Suture Passer, Frog 45° Left	801904000
EasyPass Suture Passer, Frog Crescent	801915700
EasyPass Suture Passer, Frog 70° Left	801915800
EasyPass Suture Passer, Frog 70° Right	801915900
EasyPass Suture Passer, Frog 70° Straight	RJ8100397