

rejoin 

rejoin 

Hangzhou Rejoin Mastin Medical Device Co., Ltd  
22 Xinyan Rd., Yuhang Economic Development  
Zone, Hangzhou, Zhejiang, 311100  
P. R. of China

[contact@rejoin.com](mailto:contact@rejoin.com)

[www.rejoin-medical.com](http://www.rejoin-medical.com)

 [@rejoinmedical](#)  [@rejoinmedical](#)  [@rejoinmedical](#)  [@rejoinmedical](#)


**PEEK Anchor**



## SolidPrint® Anchor, PEEK

SLAP & Bankart

- Preloaded WilSuture™ Ultra
- Size: 2.3mm
- Radiolucent and strong PEEK material
- Ergonomic handle design provides easy & quick anchor deployment



## FusionLoc® Anchor, PEEK

Rotator Cuff Repair

- Open & vented design facilitates bone ingrowth, accelerates healing
- Extended, fully-threaded anchor design provides significantly more threaded engagement, delivering greater pullout strength in poor-quality, osteoporotic bone
- Inserter reaches the tip of the anchor, spreads the stress & achieves better torsion force
- Size: 4.5mm & 5.5mm

## SolidPrint® Curved Anchor, PEEK

SLAP & Bankart & Hip Labral Repair

- Preloaded WilSuture Ultra
- Radiolucent and strong PEEK material
- Size: 2.3mm
- Curved guiding design for optimal anchor implanting angle

