

LIQUIBAND® Optima

TOPICAL SKIN ADHESIVE

liquiband.com



Take the *sting* out
of wound closure


 **LIQUIBAND®**
TOPICAL ADHESIVES

Microbial barrier
in 10 seconds¹

LIQUIBAND[®] Optima

TOPICAL SKIN ADHESIVE

LiquiBand[®] Optima Topical Skin Adhesive's unique, winged applicator allows for **effortless activation** and **controlled application**.^{3, 4}




Designed for closure of clean, fresh wounds with easily-apposed, tension-free edges



Contains a sterile blend of 90% butyl and 10% octyl cyanoacrylate that can be stored at ambient temperatures

*When applying LiquiBand[®] Optima in areas of dense hair growth, hair must be parted before applying adhesive.



**Ideal for use in any setting, including
Emergency Departments, Surgery
Centers, and Physician Offices**

- ▶ Sting-free application with **10-second dry time**^{1,13}
- ▶ Acts as a **water-resistant microbial barrier**, allowing patients to shower with physician approval immediately⁶⁻⁸
- ▶ Offers **excellent cosmesis**⁹
- ▶ Suitable for the **closure of acute wounds** including scalp, chin, facial lacerations, and punctures¹¹
- ▶ **Sloughs off naturally** within 5-14 days after application



Achieve tensile strength equal to a 4-0 suture with **just one layer** of LiquiBand® Optima¹⁰

LIQUIBAND® Optima

TOPICAL SKIN ADHESIVE



Fast, Precise, Comfortable.

Fast: close wounds in as little as 10 seconds¹

Precise: ideal applicator for facial and scalp wounds¹¹

Comfortable: no-sting adhesive for minimal patient discomfort^{12, 13}



10 SECOND
SET TIME



MICROBIAL
BARRIER



PRECISION
APPLICATOR



NO
STING



AMBIENT
STORAGE



MADE IN
THE UK

Product	Volume Per Device	Qty Per Box	Product Code
Optima	0.5g	10	OPT001
Optima Mini	0.35g	10	OPT002

Advanced Medical Solutions Limited

Langage Science Park, Western Wood Way, Plymouth PL7 5BG

www.liquiband.com | info@liquiband.com | +44 (0) 1752 209 955



Advanced Medical Solutions