



LiquiBand® Exceed™

We've drawn the line for an even, consistent wound coverage.




✓ MARKET LEADING YIELD

- ▶ 0.8ml per applicator allows for closure of wounds up to 30cm in length
- ▶ No clog technology allows intraoperative reuse for up to 90 minutes on a single patient

✓ PURE OCTYL FORMULATION

- ▶ Maximises durability and flexibility
- ▶ Provides high viscosity
- ▶ Creates intimate skin contact



✓ UNIQUE APPLICATOR TIP

- ▶ Elliptical shape allows for narrow or wide application
- ▶ Porous felt tip allows drip-free priming and even, consistent application

✓ EASY TO USE APPLICATOR



- ▶ Winged applicator allows for safe and easy activation
- ▶ Provides controlled delivery of adhesive



Better wound closure starts with a better tool


LiquiBand® Exceed™

Offers the improved attributes of precision and controlled delivery with added advantages of a new tip design.



Key Benefits

- ✓ **High Viscosity³**
 - ▶ Provides optimal flexibility
- ✓ **Controlled Expression of Adhesive⁴**
 - ▶ Applicator offers control during expression of adhesive
- ✓ **Microbial Barrier⁵**
 - ▶ Acts as a barrier against gram-positive, gram-negative and fungal microbes. Proven effectiveness against *S. aureus*, *Paeruginosa*, *E.coli*, *Candida albicans* and MRSA
- ✓ **Better Yield of Adhesive⁶**
 - ▶ 0.8 mL volume allows for greater wound coverage, up to 30cm, making it more cost effective
 - ▶ Intraoperative reuse for up to 90 minutes on a single patient
- ✓ **Maintains wound closure throughout healing process⁷**
 - ▶ The adhesive will stay on the skin for 5-10 days and sloughs off naturally



LiquiBand® Exceed™

Size	Qty Per Box	Product Code	NHS SC Code
0.8g	6	72014002	FVF022