

Contact details

If you are interested in
ELBOW PLATES
please feel free to contact us



Königsee Implantate GmbH

OT Aschau | Am Sand 4
07426 Allendorf
Germany
export@koenigsee-implantate.de
+ 49 36738 498-560

or your local sales representatives.

www.koenigsee-implantate.de

Matching screws

Cortical screw 3.5
fully threaded; self tapping;
titanium

Internal hexagon



Cortical screw 3.5
angle-stable; fully threaded;
self tapping; titanium

Internal hexagon



Cortical screw 4.5
fully threaded; self tapping;
titanium

Internal hexagon



Cortical screw 4.5
angle-stable; fully threaded;
self tapping; titanium

Internal hexagon



Distal humerus plate Y 3.5/4.5 as

Angle-stable; plate thickness head 3.2 mm; plate
thickness shaft: 4.5 mm; for screws \varnothing 3.5/4.5 mm;
8 head holes; titanium



- 5 – 8 shaft holes
- length 141 – 199 mm



ELBOW PLATES

Titanium plates range



Distal humerus Plate 3.5 as radial without extension

Angle-stable; plate thickness 3.0 mm; for screws \varnothing 3.5 mm; 4 head holes; titanium



- 4 – 10 shaft holes
- length 78 – 150 mm
- left



- 4 – 10 shaft holes
- length 78 – 150 mm
- right

Distal humerus Plate 3.5 as radial with extension

Angle-stable; plate thickness 3.0 mm; for screws \varnothing 3.5 mm; 5 head holes; titanium



- 4 – 10 shaft holes
- length 78 – 150 mm
- left



- 4 – 10 shaft holes
- length 78 – 150 mm
- right



Distal humerus plate 3.5 vas ulnar

Variable angle-stable; plate thickness head 2.9 mm; plate thickness shaft 2.6 mm; for screws \varnothing 3.5 mm; 4 head holes; titanium



- 4 – 10 shaft holes
- length 76 – 148 mm
- left



- 4 – 10 shaft holes
- length 76 – 148 mm
- right



Olecranon plate 3.5 vas with claws

Variable angle-stable; plate thickness 3.0 mm; for screws \varnothing 3.5 mm; titanium



- 8 – 14 holes
- length 76 – 151 mm



Olecranon claw plate 3.5*

Plate thickness 2.0 mm; for screws \varnothing 3.5 mm; width 10 mm; titanium



- 4 – 14 holes
- 1 long hole
- length 68 – 188 mm

* These plates are also available in steel

Olecranon claw plate 3.5 long head

Plate thickness 2.0 mm; for screws \varnothing 3.5 mm; width 11 mm; titanium



- 6 – 12 holes
- 2 – 3 long holes
- length 136 – 212 mm

Olecranon plate 3.5 vas without claws

Variable angle-stable; plate thickness 3.0 mm; for screws \varnothing 3.5 mm; titanium



- 8 – 14 holes
- length 76 – 151 mm

