

# TYMPANOMETRY ADD-ON FOR SENTIERO ADVANCED



PATH MEDICAL GmbH

Landsberger Straße 65  
82110 Germering / Germany

Tel +49 89 800 76 502 / Fax +49 89 800 76 503  
sales@pathme.de / www.pathme.de




## Tympanometry modules to „add-on“ for SENTIERO ADVANCED

The tympanometry add-on is designed to increase and complete the diagnostic inventory of your SENTIERO ADVANCED. Now you can add the powerful and proven middle ear features to your SENTIERO ADVANCED.


This unique all-in-one instrument combines OAE, Tympanometry, Pure Tone and Speech Audiometry and Evoked Potentials (ABR, ASSR, ECoChG) including both screening and diagnostic features.

From now on you can add these tympanometry options too:

Features	Screening	Screening + ETF	Screening + 1000 Hz	Diagnostic
<b>Order Number</b>	100560	100913	100689	100561
<b>Pressure Range</b>				
-300 to +300 daPa	•	•	•	•
-600 to +400 daPa				•
Manual pressure control				•
<b>Frequencies</b>				
226 Hz tone	•	•	•	•
1000 Hz tone			•	•
678 Hz, 800 Hz tone; simultaneous multi-frequency recording (all frequencies)				•
3D graph display available				•
<b>Pump Speed</b>				
50, 100, 150, 200 daPa/s; as fast as possible	•	•	•	•
<b>Extras</b>				
Children entertainment mode	•	•	•	•
Y/B/G components view (admittance, susceptance, conductance)				•
<b>Acoustic Reflex Testing</b>				
500, 1000, 2000, 3000, 4000 Hz stimulus up to 100 dB (105 dB HL for some)	•	•	•	•
Broadband, high and low-pass-Noise stimu- lus up to 90 dB HL	•	•	•	•
Ipsilateral reflex	•	•	•	•
Contralateral reflex mono insert phone mono headphone	optional #100867 #100725	optional #100867 #100725	optional #100867 #100725	•
Automatic reflex threshold	•	•	•	•
Manual reflex				•
<b>ETF - Eustachian tube function test</b>				
Intact eardrum testing		•	optional	•
Perforated eardrum testing		•	optional	•
Patulous eustachian tube testing		•	optional	•
<b>Upgradable</b>				
•	•	•	•	•
<b>Dimensions and weight</b>				
150 x 72 x 55 mm 170 grams	•	•	•	•



Portable  
Reliable  
Screening and diagnostic  
(class 2 and class 1)



Automated classification  
Multi-frequency tymp



Automated Reflex threshold  
Ipsi and contra

Additional technical specifications can be found in the detailed technical manual available online:

<https://pathme.de/support/#manuals>



Order instructions for Tympanometry Add On Hardware Type TY - MA:

# 100947-01: TY-MA incl. #100560 and a selection of ear tips

# 100947-02: TY-MA incl. #1100561 and a selection of ear tips

Please add ETF (# 100913) and / or 1000 Hz (# 100689) in addition to #100947-01 if appropriate.

CE  
0124

**Made in Germany**