

**AOHUA**



**ACD-1** |

Endoscopic CO<sub>2</sub> Insufflation Device



# Product Features

- An intelligent gas source pressure measurement device, which is designed to prevent the output flow from being too high or low.
- It can be directly connected to the carbon dioxide tubing in hospitals.
- It is compact and light to be easily installed on the endoscope cart.
- The operation is easy and convenient with a one-touch switch.

<b>Power Supply</b>	Voltage	100-240V AC
	Frequency	50/60 Hz
	Power	50 VA
<b>Specification</b>	Dimensions	327mmx180mmx228mm

 **86-21-6768 1018**

**Shanghai Aohua Photoelectricity Endoscope CO.,LTD.**  
No.66, Lane133, Guangzhong Road,Minhang District,  
Shanghai, P.R.China.(Postcode: 201108)  
Email: marketing@aohua.com  
**www.aohua.com/en**

See the Manual for details of contraindications or precautions