





Product Features

- It is easy to operate, can provide powerful flushes by using a foot switch.
- The flushing system is completely closed and is directly connected to a bottle of sterile water, which makes it convenient and fast to be used.
- During the diagnosis with an endoscope, irrigation should be carried out at any time to wash the digestive tract mucosa and dissection site and to clean the surgical field of vision.
- Automatically stops working 20 seconds after starting up: AFP-1 will stop working after pumping water for 20 seconds.
- The maximum flow rate can be up to 450 ml per minute and the flow rate is displayed in sections.
- Compact and light: it is light and compact and can be directly placed on a cart.
- The pump head is removable: the design of removable pump head facilitates the assembly.

	Voltage	100-240V AC
Power Supply	Frequency	50/60 Hz
	Power	50 VA
Specification	Continuous Water Supply	20s



Shanghai Aohua Photoelectricity Endoscope CO.,LTD.
No.66, Lane133, Guangzhong Road,Minhang District,
Shanghai, P.R.China.(Postcode: 201108)
Email: marketing@aohua.com

www.aohua.com/en