



CRUISER BRACE



The Rhino Cruiser is a hip abduction brace fabricated from polypropylene, for infants to 3 years of age with hip dysplasia. The Rhino hold the hips in abduction, whilst allowing the child to walk. Can be used full time or during the night only. Often the Rhino is used after treatment with a pavlik harness, once the child begins to walk.

Sizes	Waist	Item Number
Preemie	12" - 14"	27090
Ex-Small	14" - 16"	27091
Small	16" - 18"	27092
Medium	18" - 20"	27093
Large	20" - 22"	27094
Ex-Large	22" - 24"	27095

- Durable plastic shell lined with closed-cell foam
- Lightweight and open
- Easy to apply
- Easy to clean – foam will not contain body fluids or excrement
- 90 degree hip flexion, 60 degree abduction

