




UroDapter

UroDapter: Revolutionizing Intravesical Instillation



The UroDapter is a groundbreaking catheter-free syringe adapter transforming the approach to intravesical instillations. Developed by UroSystem, our mission is to enhance the quality of life for individuals suffering from bladder diseases by introducing practice-changing medical solutions.

Uniquely designed for use in both male and female patients, the UroDapter ensures a simpler, safer, and more comfortable solution to a wide array of urological conditions, including but not limited to:

- Interstitial Cystitis / Bladder Pain Syndrome
- Recurrent Urinary Tract Infections
- Bladder Cancer
- Post-Irradiation Cystitis
- Post-chemotherapy Cystitis

Catheter-Free Convenience



Designed for both male and female patients, the UroDapter offers a non-invasive alternative to catheters, optimizing efficacy while minimizing the risk of complications.

Enhanced Patient Comfort



By eliminating the need for catheterization, the UroDapter reduces patient discomfort, anxiety, and the psychological distress associated with catheter-based instillations.

Superior Safety Profile



Reduces the risk of urinary tract infections (UTIs), urethral bleeding and strictures, ensuring a safer treatment option for intravesical therapy.

Efficiency in Treatment



Enables quick and easy administration, saving valuable time for healthcare professionals while ensuring effective delivery of therapeutic agents for bladder and urethral treatment.

Technical Specifications



- **Ergonomic Design:** Features a short, rounded tip for minimal urethral insertion, an isolating collar for leak-proof instillation, and a ribbed handling grip for ease of use.
- **Compatibility:** Designed to connect seamlessly with both Luer slip and Luer lock syringes, facilitating the use of a wide range of medicinal solutions.

For detailed step-by-step user manuals please visit <https://app.urodapter.com/>.

Implementing UroDapter in Your Practice

Incorporating UroDapter into your clinical practice is straightforward. Our support includes comprehensive training materials, instructional videos, and a dedicated customer service team to assist with any inquiries. Elevate your standard of care with UroDapter, and join us in setting new benchmarks in urological treatment.

Contact Us

For more information please contact your local distributor or contact us directly.

UroSystem:
Elevating Urological Care

